## SISTEMAS DE IMPLANTES IMPLANT SYSTEMS

[www.gmidental.com](http://www.gmidental.com)

avantgard

avantgard  
**PEAK**

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### KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

### IMPLANT SYSTEMS

- **gmi® avantgard** dental implants 8
- Insertion sequence for **gmi® avantgard** implants 9
- **gmi® avantgard** implants range 10
- Drilling sequences for **gmi® avantgard** implants 12
- **gmi® avantgard PEAK** dental implants 14
- Insertion sequence for **gmi® avantgard PEAK** implants 15
- **gmi® avantgard PEAK** implants range 16
- Drilling sequences for **gmi® avantgard PEAK** implants 17

### GMI® AVANTGARD ATTACHMENTS

- Healing abutments 19
- Implant carrier abutments 20
- Temporary cylinders 21
- Angled abutments 22
- Straight abutments 23
- Titanium base 24
- Multi-esthetic abutments 26
- Bioesthetic bases 29
- CrCo pre-milled blanks 30
- Replicas 31
- Impression coping 32
- Clinic and laboratory screws 33
- MEA and BE screws 34
- Recommended tightening torques 35

### GMI® SURGICAL KITS

- **gmi® frontier/avantgard** surgical kit 37
- **gmi® avantgard** guided surgery kit 38
- **gmi® PEAK** surgical kit 39
- **gmi® avantgard PEAK** guided surgery kit 40

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- The **CE 2797** marking applies to medical devices certified through conformity assessment by a Notified Body under Regulation (EU) 2017/745.

- The **CE** marking (without an identification number) applies to self-certified medical devices under Regulation (EU) 2017/745.




- Neither of these symbols applies to products that fall outside the scope of Regulation (EU) 2017/745 for medical devices.

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To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.





SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

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The **gmi**® **avantgard** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

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The **gmi**® **avantgard** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi**® recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi**® **avantgard** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi**® recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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

To prevent prior deformation of the clinic screw, **gmi**® recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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
In the chart are detailed all the symbols appearing on the gmi® avantgard implant system's labelling and packaging and their corresponding description.

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\*Minimum requirements for downloading and opening the electronic instructions for use:

- Stable internet connectivity.
- A device with a modern operating system that supports PDF viewing.
- Adequate storage space on your device.
- A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.
- A PDF reader application installed on your device.

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
gmi® avantgard conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high primary stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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⑥ Rosca interior  
NP-RP: M1,6 x 0,35 mm  
WP: M2,0 x 0,40 mm

⑦ Cuerpo del implante  
Paso: 1,60 mm

① Platforms  
NP: narrow  
RP: standard  
WP: wide

② Length

③ Self threading millings  
Variable length  
3 millings at 120°  
Atraumatic tip

④ Cortical thread

⑤ Internal grooves  
6 x 60°

⑥ Internal thread  
NP-RP: M1.6 x 1.35 mm  
WP: M2.0 x 0.40 mm

⑦ Implant body  
Pitch: 1.60 mm

|   | NP                          | RP                        | RP                            | WP                            | WP                                |
|---|-----------------------------|---------------------------|-------------------------------|-------------------------------|-----------------------------------|
| ① | NP                          | RP                        | RP                            | WP                            | WP                                |
| ② | De 10 a 15<br>From 10 to 15 | De 8 a 15<br>From 8 to 15 | De 6,5 a 13<br>From 6.5 to 13 | De 6,5 a 13<br>From 6.5 to 13 | De 6,5 a 11,5<br>From 6.5 to 11.5 |
| ④ | ø3,30                       | ø4,00                     | ø4,50                         | ø5,00                         | ø6,00                             |
| ⑦ | ø3,30                       | ø3,75                     | ø4,25                         | ø4,75                         | ø5,75                             |

## Protocolo de inserción de implantes gmi® avantgard - Insertion sequence for gmi® avantgard implants

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gmi® avantgard conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the gmi® avantgard implants consists of a titanium support which houses the dental implant, prepared to be inserted directly into the mouth, using the implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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\* Consult the instructions for use in electronic format.


1. Check sticker indicating implant diameter and length and open the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the insertion wrench facing the grooves of the wrench with those of the implant.
7. Check that the wrench is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the subcrestal level. Remove the wrench by making an axial movement.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove wrench and suture the incision.

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

The following is a description of the drilling sequences for the different conical connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped drill: 1200 - 1500 rpm.
- Pilot drill: 700 - 900 rpm.
- Step drills:
  - Ø3.00 & Ø3.50 mm → 400 - 700 rpm.
  - Ø4.00 & Ø5.10 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.

For drilling, external cooling with saline solution is required.


| avantgard Ø3,30           |                              |                        |           |
|---------------------------|------------------------------|------------------------|-----------|
| Descripción - Description |                              |                        | Ref.      |
| 1                         | Fresa lanza                  | Lance-shaped drill     | KYFoC1225 |
| 2                         | Fresa piloto Ø2,50           | Pilot drill Ø2.50      | KYFoC2250 |
| 3                         | Fresa cortical Ø3,20         | Cortical drill Ø3,20   | KYFoC4332 |
| 4                         | Fresa escalonada Ø2,50/Ø3,00 | Step drill Ø2.50/Ø3.00 | KYFoC5130 |


| avantgard Ø3,75           |                              |                        |           |
|---------------------------|------------------------------|------------------------|-----------|
| Descripción - Description |                              |                        | Ref.      |
| 1                         | Fresa lanza                  | Lance-shaped drill     | KYFoC1225 |
| 2                         | Fresa piloto Ø2,50           | Pilot drill Ø2.50      | KYFoC2250 |
| 3                         | Fresa escalonada Ø2,50/Ø3,00 | Step drill Ø2.50/Ø3.00 | KYFoC5130 |
| 4                         | Fresa cortical Ø4,00         | Cortical drill Ø4.00   | KYFoC4340 |
| 5                         | Fresa escalonada Ø3,00/Ø3,50 | Step drill Ø3.00/Ø3.50 | KYFoC5135 |

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

gmi® avantgard PEAK conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the conical core and progressive thread profile facilitate implant insertion and act as a bone compactor allowing high primary stability to be obtained in sparse bones.

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② Cortical thread  
1.50 mm

③ Length

④ Self threading millings  
Variable length  
2 millings at 180°  
Atraumatic tip  
Cut in both directions

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gmi® avantgard PEAK conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the gmi® avantgard PEAK implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single Avantgard NP-RP implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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\* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the NP-RP insertion wrench facing the grooves of the wrench with those of the implant.
7. Check that the wrench is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the subcrestal level. Remove the wrench by making an axial movement.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove wrench and suture the incision.

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

The following is a description of the drilling sequences for the different conical connection implant models, as well as the recommended drilling conditions for their use:

- Initial drill: 1000 rpm.
- Pilot drill: 700 - 900 rpm.
- Step and conical drills:
  - ∅2.50 / ∅2.80 mm → 500 - 700 rpm.
  - ∅3.20 / ∅3.65 mm → 500 - 700 rpm.
  - ∅3.70 / ∅4.15 mm → 400 - 700 rpm.
  - ∅4.15 / ∅4.65 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.
- For drilling, external cooling with saline solution is required.


| avantgard PEAK Ø3,30      |                                 |                           |           |
|---------------------------|---------------------------------|---------------------------|-----------|
| Descripción - Description |                                 | Ref.                      |           |
| 1                         | Fresa inicial                   | Initial drill             | KYFOC1401 |
| 2                         | Fresa piloto Ø2,50 mm           | Pilot drill Ø2.50         | KYFOC2250 |
| 3                         | Fresa escalonada Ø2,50/Ø2,80 mm | Step drill Ø2.50/Ø2.80 mm | KYFOC5129 |
| 4                         | Fresa cortical Ø3,20 mm         | Cortical drill Ø3.20      | KYFOC4333 |


| avantgard PEAK Ø3,75      |                                 |                           |           |
|---------------------------|---------------------------------|---------------------------|-----------|
| Descripción - Description |                                 | Ref.                      |           |
| 1                         | Fresa inicial                   | Initial drill             | KYFOC1401 |
| 2                         | Fresa piloto Ø2,50 mm           | Pilot drill Ø2.50         | KYFOC2250 |
| 3                         | Fresa escalonada Ø2,50/Ø2,80 mm | Step drill Ø2.00/Ø2.80 mm | KYFOC5129 |
| 4                         | Fresa cónica Ø3,20 mm           | Conical drill Ø3.20 mm    | KYFOC7032 |
| 5                         | Fresa cortical Ø3,65 mm         | Cortical drill Ø3.65      | KYFOC4337 |
| 6                         | Fresa escalonada Ø3,20/Ø3,65 mm | Step drill Ø3.20/Ø3.65 mm | KYFOC5137 |

## Secuencias fresado para implantes gmi® avantgard PEAK - Drilling sequences for gmi® avantgard PEAK implants

| avantgard PEAK $\varnothing 4,25$ |   |   |
|-----------------------------------|---|---|
|                                   | Descripción - Description                               | Ref.  |
| 1                                 | Fresa inicial   | Initial drill                                     |
| 2                                 | Fresa piloto $\varnothing 2,50$ mm                      | Pilot drill $\varnothing 2.50$                    |
| 3                                 | Fresa escalonada $\varnothing 2,50/\varnothing 2,80$ mm | Step drill $\varnothing 2.00/\varnothing 2.80$ mm |
| 4                                 | Fresa cónica $\varnothing 3,20$ mm                      | Conical drill $\varnothing 3.20$ mm               |
| 5                                 | Fresa cónica $\varnothing 3,70$ mm                      | Conical drill $\varnothing 3.70$ mm               |
| 6                                 | Fresa cortical $\varnothing 4,15$ mm                    | Cortical drill $\varnothing 4.15$ mm              |
| 7                                 | Fresa escalonada $\varnothing 3,70/\varnothing 4,15$ mm | Step drill $\varnothing 3.70/\varnothing 4.15$ mm |

| avantgard PEAK $\varnothing 4,75$ |   |   |
|-----------------------------------|---|---|
|                                   | Descripción - Description                               | Ref.  |
| 1                                 | Fresa inicial   | Initial drill                                     |
| 2                                 | Fresa piloto $\varnothing 2,50$ mm                      | Pilot drill $\varnothing 2.50$                    |
| 3                                 | Fresa escalonada $\varnothing 2,50/\varnothing 2,80$ mm | Step drill $\varnothing 2.00/\varnothing 2.80$ mm |
| 4                                 | Fresa cónica $\varnothing 3,20$ mm                      | Conical drill $\varnothing 3.20$ mm               |
| 5                                 | Fresa cónica $\varnothing 3,70$ mm                      | Conical drill $\varnothing 3.70$ mm               |
| 6                                 | Fresa cónica $\varnothing 4,20$ mm                      | Conical drill $\varnothing 4.20$ mm               |
| 7                                 | Fresa cortical $\varnothing 4,65$ mm                    | Cortical drill $\varnothing 4.65$ mm              |
| 8                                 | Fresa escalonada $\varnothing 4,15/\varnothing 4,65$ mm | Step drill $\varnothing 4.15/\varnothing 4.65$ mm |

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Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Available in emergency heights from 1.50 to 4.50 mm.
- Apical bit to facilitate the start of threading.
- Manual tightening torque: Max. 15 N-cm.

| Plataforma<br>Platform | Ø cabeza<br>Ø head | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|--------------------|----------------------------|-----------|
| NP ●●                  | Ø 4,00             | 1,50                       | KDPoC4012 |
|                        |                    | 2,50                       | KDPoC4013 |
|                        |                    | 3,50                       | KDPoC4014 |
|                        |                    | 4,50                       | KDPoC4015 |
| RP ●●                  | Ø 5,00             | 1,50                       | KDPoC4016 |
|                        |                    | 2,50                       | KDPoC4017 |
|                        |                    | 3,50                       | KDPoC4018 |
|                        |                    | 4,50                       | KDPoC4019 |

| Plataforma<br>Platform | Ø cabeza<br>Ø head | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|--------------------|----------------------------|-----------|
| NP ●●                  | Ø 5,00             | 1,50                       | KDPoC4001 |
|                        |                    | 2,50                       | KDPoC4002 |
|                        |                    | 3,50                       | KDPoC4003 |
|                        |                    | 4,50                       | KDPoC4004 |
| RP ●●                  | Ø 6,00             | 2,50                       | KDPoC4009 |
|                        |                    | 3,50                       | KDPoC4010 |
|                        |                    | 4,50                       | KDPoC4011 |

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
When the use of a custom abutment is required, implant carrier abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, implant carrier abutments can be used like closed tray impression coping.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 3.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping.
- Premounted clinic screw included.
- Screw tightening torque: 25 N-cm.

| Plataforma<br>Platform | Ø hombro<br>Ø shoulder | Ref.      |
|------------------------|------------------------|-----------|
| RP ●●                  | Ø 3,75 - Ø 4,25        | KDBoF4001 |
| WP ●                   | Ø 4,75                 | KDBoF4002 |
|                        | Ø 5,75                 | KDBoF4003 |

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
When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary abutment is used, defined as attachment that once attached to the crown of the implant with clinic screw, they will allow a temporary restoration to be carried out while the definitive restoration is being manufactured.

### Key features

- Machined in high strength technical plastic.
- NP abutments machined in Grade V titanium as per ISO-5832 standard
- Rotary version only for multiple restorations.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw included.
- Manual tightening torque: Max. 15 N·cm.

| Plataforma<br>Platform | Tipo conexión<br>Connection type | Ref.      |
|------------------------|----------------------------------|-----------|
| RP ● ●                 | Anti-rotatoria • Non-rotary      | KDHOF4001 |
|                        | Rotatoria • Rotary               | KDHOF4002 |
| WP ●                   | Anti-rotatoria • Non-rotary      | KDHOF4003 |
|                        | Rotatoria • Rotary               | KDHOF4004 |

| Plataforma<br>Platform | Altura (mm)<br>Height (mm) | Tipo conexión<br>Connection type | Ref.      |
|------------------------|----------------------------|----------------------------------|-----------|
| NP ● ●                 | 1,50                       | Anti-rotatoria • Non-rotary      | KDHOF4009 |
|                        | 3,00                       | Anti-rotatoria • Non-rotary      | KDHOF4010 |
|                        | 1,50                       | Rotatoria • Rotary               | KDHOF4011 |
|                        | 3,00                       | Rotatoria • Rotary               | KDHOF4012 |

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw included.
- Abutment shoulder adapted to the shape of the gum.
- Screw tightening torque: 25 N-cm.

| Plataforma<br>Platform | Ang. (°) | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|----------|----------------------------|-----------|
| NP ●●                  | 10°      | 2,50                       | KDMoF4039 |
|                        | 10°      | 4,50                       | KDMoF4040 |
|                        | 20°      | 2,50                       | KDMoF4041 |
|                        | 20°      | 4,50                       | KDMoF4042 |

| Plataforma<br>Platform | Ang. (°) | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|----------|----------------------------|-----------|
| RP ●●                  | 15°      | 2,50                       | KDMoF4013 |
|                        | 15°      | 4,50                       | KDMoF4014 |
|                        | 20°      | 2,50                       | KDMoF4015 |
|                        | 20°      | 4,50                       | KDMoF4016 |

| Plataforma<br>Platform | Ang. (°) | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|----------|----------------------------|-----------|
| WP ●                   | 15°      | 2,50                       | KDMoF4017 |
|                        | 15°      | 4,50                       | KDMoF4018 |
|                        | 20°      | 2,50                       | KDMoF4019 |
|                        | 20°      | 4,50                       | KDMoF4020 |

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
Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

### Key features


- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Available only in non-rotary version.
- Clinic screw included.
- Available in four emergency heights from 1.50 to 4.50 mm.
- Screw tightening torque: 25 N-cm.

| Plataforma<br>Platform | Ø hombro<br>Ø shoulder | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|------------------------|----------------------------|-----------|
| NP ●●                  | Ø 4,00                 | 1,50                       | KDMoF4035 |
|                        |                        | 2,50                       | KDMoF4036 |
|                        |                        | 3,50                       | KDMoF4037 |
|                        |                        | 4,50                       | KDMoF4038 |

| Plataforma<br>Platform | Ø hombro<br>Ø shoulder | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|------------------------|----------------------------|-----------|
| RP ●●                  | Ø 4,00                 | 1,50                       | KDMoF4001 |
|                        |                        | 2,50                       | KDMoF4002 |
|                        |                        | 3,50                       | KDMoF4003 |
|                        |                        | 4,50                       | KDMoF4004 |
|                        | Ø 5,00                 | 1,50                       | KDMoF4005 |
|                        |                        | 2,50                       | KDMoF4006 |
|                        |                        | 3,50                       | KDMoF4007 |
|                        |                        | 4,50                       | KDMoF4008 |

| Plataforma<br>Platform | Ø hombro<br>Ø shoulder | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|------------------------|----------------------------|-----------|
| WP ●                   | Ø 6,00                 | 1,50                       | KDMoF4009 |
|                        |                        | 2,50                       | KDMoF4010 |
|                        |                        | 3,50                       | KDMoF4011 |
|                        |                        | 4,50                       | KDMoF4012 |

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Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw included.
- Available in emergency height of 1.50 mm and 3.00 mm.
- Non-rotary titanium base with indexation.
- Screw tightening torque: 25 N-cm.
- Scan directly on Titanium base or alternative scan on Scan abutment.

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

Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.50 to 4.50 mm.
- One-piece straight abutments in rotary version.
- Tightened with hex wrench: B/F 2.00 mm (KYLOCO149).
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.
- Tightening torques: ME straight abut. (30 N-cm), ME angled abut. (20 N-cm) and healing caps (Manual: Max. 15 N-cm).

### Pilares multi-estético rectos Multi-esthetic straight abutments

| Plataforma<br>Platform | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|----------------------------|-----------|
| NP ●●                  | 1,50                       | KDGoF4019 |
|                        | 2,50                       | KDGoF4020 |
|                        | 3,50                       | KDGoF4021 |
|                        | 4,50                       | KDGoF4022 |
| RP ●●                  | 1,50                       | KDGoF4001 |
|                        | 2,50                       | KDGoF4002 |
|                        | 3,50                       | KDGoF4003 |
|                        | 4,50                       | KDGoF4004 |
| WP ●                   | 1,50                       | KDGoF4005 |
|                        | 2,50                       | KDGoF4006 |
|                        | 3,50                       | KDGoF4007 |
|                        | 4,50                       | KDGoF4008 |

## Pilares multi-estético gmi® avantgard - gmi® avantgard multi-esthetic abutments

### Pilares multi-estético inclinados Multi-esthetic angled abutments

| Plataforma<br>Platform | Ang. (°) | Altura (mm)<br>Height (mm) | Ref.      |
|------------------------|----------|----------------------------|-----------|
|                        | 17°      | 2,50                       | KDGoF4009 |
| RP ●●                  | 17°      | 3,50                       | KDGoF4010 |
|                        | 30°      | 4,00                       | KDGoF4011 |


### Pilares cicatrización PME Healing abutments MEA

| Plataforma<br>Platform | Ø Diámetro (mm)<br>Ø Diameter (mm) | Ref.      |
|------------------------|------------------------------------|-----------|
| NP - RP - WP ●●        | Ø 4,80                             | KDPoF0035 |
|                        | Ø 6,00                             | KDPoF0037 |

### Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

| Plataforma<br>Platform | Ref.      |
|------------------------|-----------|
| NP - RP - WP ●●        | KDSoF0015 |

### Tomas impresión cubeta cerrada (CC) PME Closed tray (CT) impression coping MEA

| Plataforma<br>Platform | Cantidad (Ud.)<br>Quantity (Un.) | Ref.      |
|------------------------|----------------------------------|-----------|
| NP - RP - WP ●●        | 1                                | KDSoF0017 |
|                        | 5                                | KDSoF5001 |

## Pilares multi-estético gmi® avantgard - gmi® avantgard multi-esthetic abutments

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Once the insertion of the implant in the bone bed has been completed, the bioesthetic base is used, which is defined as the machined attachment that, once fixed to the coronal part of the implant, allows single restorations at the abutment level.

Placing the bioesthetic base during surgery allows the soft tissues to heal around it.

Once the tissue has healed, it serves as the basis for the rehabilitation, avoiding the connection and disconnection of the components to the implant and, therefore, keeping the biological seal intact.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 2.50 to 4.50 mm.
- Straight two-piece bases, including the screw.
- Tightened with hex wrench: B/F 2.00 mm (KYLOCo149).
- Impression coping with screw included.
- Castable.
- Tightening torques: BE base (25 N-cm) and healing abutment (Manual: Max. 15 N-cm).

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
CrCo pre-milled blanks consist of a block with the implant connection machined, which allows an excellent fit between the implant and the abutment.

They are indicated to produce custom one-piece abutments using CAD/CAM allowing direct application of ceramic.

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

Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.


| Réplicas de implante<br>Implant replicas |           |
|--|-----------|
| Plataforma<br>Platform                   | Ref.      |
| NP ●●                                    | KDRoC4005 |
| RP ●●                                    | KDRoC4001 |
| WP ●                                     | KDRoC4002 |

| Réplicas digitales<br>Digital replicas |           |
|--|-----------|
| Plataforma<br>Platform                 | Ref.      |
| NP ●●                                  | KDRoF4006 |
| RP ●●                                  | KDRoF4003 |
| WP ●                                   | KDRoF4004 |

| Réplicas para pilares multi-estético<br>MEA replicas |           |
|--|-----------|
| Plataforma<br>Platform                               | Ref.      |
| NP - RP - WP ●●                                      | KDRoC0012 |

| Réplicas digitales para pilares multi-estético<br>MEA digital replicas |           |
|--|-----------|
| Plataforma<br>Platform   | Ref.      |
| NP - RP - WP ●●  | KDRoF0014 |

## Réplicas gmi® avantgard - gmi® avantgard replicas

| Réplica digital BE<br>BE digital replica |           |
|--|-----------|
| Plataforma<br>Platform                   | Ref.      |
| RP - WP ●●                               | KDRoF0016 |

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
Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, is used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

### Key features


- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- Manual tightening torque: Max. 15 N-cm.

| Tomas impresión cubeta abierta (CA)<br>Open tray (OT) impression coping |                   |           |
|---|-------------------|-----------|
| Plataforma<br>Platform  | Tornillo<br>Screw | Ref.      |
| NP ●●   | Largo • Long      | KDS0F4005 |
| RP ●●   | Largo • Long      | KDS0F4001 |
|   | Corto • Short     | KDS0F4003 |
| WP ●  | Largo • Long      | KDS0F4002 |
|   | Corto • Short     | KDS0F4004 |

## Tomas de impresión gmi® avantgard - gmi® avantgard impression coping

### Tomas impresión cubeta cerrada (CC) plástico transportador Closed tray (CT) plastic impression coping implant-carrier

| Plataforma<br>Platform | Ref.           |
|------------------------|----------------|
| NP ●●                  | KDSOF4007      |
|                        | KDSOF0019 (x5) |
| RP - WP ●●             | KDSOF0014 (x5) |

## Tornillos clínica y laboratorio gmi® avantgard - gmi® avantgard clinic and laboratory screws

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
Clinic screws are defined as machined devices used for the final anchoring of the abutments or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of abutments to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

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### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

| Tornillo clínica<br>Clinic screw         |                                     |           |
|--|-------------------------------------|-----------|
| Plataforma<br>Platform                   | Hex. llave (mm)<br>Hex. wrench (mm) | Ref.      |
| NP - RP ●●                               | ⊕ 1,20 mm                           | KDToC4003 |
|  | ⊕ TX                                | KDToC4021 |
| WP ●                                     | ⊕ 1,20 mm                           | KDToC4004 |
|  | ⊕ TX                                | KDToC4022 |
| Tornillo laboratorio<br>Laboratory screw |                                     |           |
| Plataforma<br>Platform                   | Hex. llave (mm)<br>Hex. wrench (mm) | Ref.      |
| NP - RP ●●                               | ⊕ 1,20 mm                           | KDToC4005 |
|  | ⊕ TX                                | KDToC4023 |
| WP ●                                     | ⊕ 1,20 mm                           | KDToC4006 |
|  | ⊕ TX                                | KDToC4024 |

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MEA and BE clinical screws are defined as machined devices used for the final anchoring of the ME and BE abutments, while laboratory screws are the machined devices used for the temporary anchoring of ME and BE abutments to the replica during the fabrication of the prostheses in the laboratory.

### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm and 2.00 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

| Tornillo PME<br>MEA screw |                                     |           |
|---------------------------|-------------------------------------|-----------|
| Plataforma<br>Platform    | Hex. llave (mm)<br>Hex. wrench (mm) | Ref.      |
| NP - RP - WP ● ●          | ○ 1,20 mm                           | KDToC0056 |
|                           | ⊛ TX                                | KDToC0070 |

| Tornillo laboratorio PME<br>MEA laboratory screw |                                     |           |
|--|-------------------------------------|-----------|
| Plataforma<br>Platform                           | Hex. llave (mm)<br>Hex. wrench (mm) | Ref.      |
| NP - RP - WP ● ●                                 | ○ 1,20 mm                           | KDToC0057 |
|  | ⊛ TX                                | KDToC0072 |

| Tornillo BE<br>BE screw |                                     |           |
|-------------------------|-------------------------------------|-----------|
| Plataforma<br>Platform  | Hex. llave (mm)<br>Hex. wrench (mm) | Ref.      |
| RP - WP ● ●             | ○ 1,20 mm                           | KDToC4016 |
|                         | ⊛ TX                                | KDToC0074 |

| Tornillo laboratorio BE<br>BE laboratory screw |                                     |           |
|--|-------------------------------------|-----------|
| Plataforma<br>Platform                         | Hex. llave (mm)<br>Hex. wrench (mm) | Ref.      |
| RP - WP ● ●                                    | ○ 1,20 mm                           | KDToC4020 |
|  | ⊛ TX                                | KDToC0075 |

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


gmi® frontier/avantgard (Ref. KYCoF3053) surgical kit consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of gmi frontier and avantgard dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi® surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134°C.

To complement this kit, gmi® has the PEAK drill Kit (Ref. KYCoF3056 - *see instrumental catalogue*), allowing it to be easily expanded for the placement of 2 additional implant systems (frontier PEAK and avantgard PEAK).

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Ref. KYCoF3059 - Kit for gmi® avantgard system only, without the drill and the wrenches exclusive to the frontier system (*see table, highlighted in red*).

| Descripción - Description                       |  | Un.      | Ref.             |
|---|--|----------|------------------|
| Fresa lanza                                     | Lance-shaped drill                                 | 1        | KYFoC1225        |
| Fresa piloto ø2,50 mm                           | Pilot drill ø2.50 mm                               | 1        | KYFoC2250        |
| Fresa escalonada ø2,50 / ø3,00 mm               | Step drill ø2.50 / ø3.00 mm                        | 1        | KYFoC5130        |
| Fresa escalonada ø3,00 / ø3,50 mm               | Step drill ø3.00 / ø3.50 mm                        | 1        | KYFoC5135        |
| Fresa escalonada ø3,50 / ø4,00 mm               | Step drill ø3.50 / ø4.00 mm                        | 1        | KYFoC5140        |
| Fresa escalonada ø4,00 / ø4,50 mm               | Step drill ø4.00 / ø4.50 mm                        | 1        | KYFoC5145        |
| <b>Fresa cortical ø3,65 mm</b>                  | <b>Cortical drill ø3.65 mm</b>                     | <b>1</b> | <b>KYFoC4336</b> |
| Fresa cortical ø3,20 mm                         | Cortical drill ø3.20 mm                            | 1        | KYFoC4332        |
| Fresa cortical ø4,00 mm                         | Cortical drill ø4.00 mm                            | 1        | KYFoC4340        |
| Fresa cortical ø4,50 mm                         | Cortical drill ø4.50 mm                            | 1        | KYFoC4345        |
| Fresa cortical ø4,90 mm                         | Cortical drill ø4.90 mm                            | 1        | KYFoC4349        |
| Fresa cortical ø5,90 mm                         | Cortical drill ø5.90 mm                            | 1        | KYFoC4359        |
| Llave carraca IP                                | TI ratchet wrench                                  | 1        | KYLoF0113        |
| <b>Llave-transportador HEX-2,45 mm IP corta</b> | <b>Short HEX-2.45 mm implant carrier TI wrench</b> | <b>1</b> | <b>KYLoF0096</b> |
| <b>Llave-transportador HEX-2,45 mm IP larga</b> | <b>Long HEX-2.45 mm implant carrier TI wrench</b>  | <b>1</b> | <b>KYLoF0111</b> |
| <b>Llave-transportador HEX-2,45 mm CA</b>       | <b>HEX-2.45 mm implant carrier HP wrench</b>       | <b>1</b> | <b>KYLoF0091</b> |
| Llave Avantgard NP/RP IP corta                  | Short Avantgard NP/RP TI wrench                    | 1        | KYLoF0157        |
| Llave Avantgard NP/RP IP larga                  | Long Avantgard NP/RP TI wrench                     | 1        | KYLoF0162        |
| Llave Avantgard WP IP corta                     | Short Avantgard WP TI wrench                       | 1        | KYLoF0168        |
| Llave Avantgard WP IP larga                     | Long Avantgard WP TI wrench                        | 1        | KYLoF0169        |
| Llave Avantgard NP/RP CA corta                  | Short Avantgard NP/RP HP wrench                    | 1        | KYLoF0151        |
| Llave Avantgard WP CA corta                     | Short Avantgard WP HP wrench                       | 1        | KYLoF0166        |
| Llave carraca corta HEX-1,20 mm                 | Short HEX-1.20 mm ratchet wrench                   | 1        | KYLoF0128        |
| Llave carraca larga HEX-1,20 mm                 | Long HEX-1.20 mm ratchet wrench                    | 1        | KYLoF0129        |
| Alargador fresas                                | Drill extender                                     | 1        | KYLoC0023        |
| Paralelizador ø2,00 / ø2,80 mm                  | Paralleling pin ø2.00 / ø2.80 mm                   | 3        | KYLoC0078        |
| Indicador profundidad ø2,50 mm                  | Depth gauge ø2.50 mm                               | 1        | KYLoC0165        |

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gmi® avantgard guided surgery kit (Ref. KYCoF3049) consists of an autoclavable case made of high strength resin, containing all the necessary material for performing totally guided surgeries with gmi® avantgard implants of  $\varnothing 3.30$  /  $\varnothing 3.75$  /  $\varnothing 4.25$  and  $\varnothing 4.75$  mm and lengths between 8 and 13 mm.

To complement this kit, gmi® has the PEAK guided surgery drill Kit (Ref. KYCoF3052 - *see instrumental catalogue*), allowing it to be easily expanded for the placement of 1 additional implant system (avantgard PEAK).

| Descripción - Description                                     |   | Un. | Ref.      |
|---|---|-----|-----------|
| Fresa Pin $\varnothing 1,50$ mm                               | Pin drill $\varnothing 1.50$ mm                         | 1   | KYFoC6060 |
| Trepante gingival CA $\varnothing 5,10$ mm                    | HP gingival punch $\varnothing 5.10$ mm                 | 1   | KYFoC6066 |
| Fresa inicial $\varnothing 2,5$ mm                            | Initial drill $\varnothing 2.5$ mm                      | 1   | KYFoC6101 |
| Fresa piloto DLC $\varnothing 2,5$ L8 mm                      | DLC pilot drill $\varnothing 2.5$ L8 mm                 | 1   | KYFoC6201 |
| Fresa piloto DLC $\varnothing 2,5$ L10 mm                     | DLC pilot drill $\varnothing 2.5$ L10 mm                | 1   | KYFoC6202 |
| Fresa piloto DLC $\varnothing 2,5$ L11,5 mm                   | DLC pilot drill $\varnothing 2.5$ L11.5 mm              | 1   | KYFoC6203 |
| Fresa piloto DLC $\varnothing 2,5$ L13, mm                    | DLC pilot drill $\varnothing 2.5$ L13 mm                | 1   | KYFoC6204 |
| Fresa piloto DLC $\varnothing 2,5$ L15 mm                     | DLC pilot drill $\varnothing 2.5$ L15 mm                | 1   | KYFoC6205 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L8 mm    | DLC step drill $\varnothing 2.5/\varnothing 3$ L8 mm    | 1   | KYFoC6506 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L10 mm   | DLC step drill $\varnothing 2.5/\varnothing 3$ L10 mm   | 1   | KYFoC6507 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L11,5 mm | DLC step drill $\varnothing 2.5/\varnothing 3$ L11.5 mm | 1   | KYFoC6508 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L13 mm   | DLC step drill $\varnothing 2.5/\varnothing 3$ L13 mm   | 1   | KYFoC6509 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L15 mm   | DLC step drill $\varnothing 2.5/\varnothing 3$ L15 mm   | 1   | KYFoC6510 |
| Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L8 mm    | DLC step drill $\varnothing 3/\varnothing 3.5$ L8 mm    | 1   | KYFoC6511 |
| Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L10 mm   | DLC step drill $\varnothing 3/\varnothing 3.5$ L10 mm   | 1   | KYFoC6512 |
| Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L11,5 mm | DLC step drill $\varnothing 3/\varnothing 3.5$ L11.5 mm | 1   | KYFoC6513 |
| Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L13 mm   | DLC step drill $\varnothing 3/\varnothing 3.5$ L13 mm   | 1   | KYFoC6514 |
| Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L15 mm   | DLC step drill $\varnothing 3/\varnothing 3.5$ L15 mm   | 1   | KYFoC6515 |
| Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L8 mm    | DLC step drill $\varnothing 3.5/\varnothing 4$ L8 mm    | 1   | KYFoC6516 |
| Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L10 mm   | DLC step drill $\varnothing 3.5/\varnothing 4$ L10 mm   | 1   | KYFoC6517 |
| Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L11,5 mm | DLC step drill $\varnothing 3.5/\varnothing 4$ L11.5 mm | 1   | KYFoC6518 |
| Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L13 mm   | DLC step drill $\varnothing 3.5/\varnothing 4$ L13 mm   | 1   | KYFoC6519 |
| Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L15 mm   | DLC step drill $\varnothing 3.5/\varnothing 4$ L15 mm   | 1   | KYFoC6520 |
| Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L8 mm    | DLC step drill $\varnothing 4/\varnothing 4.5$ L8 mm    | 1   | KYFoC6521 |
| Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L10 mm   | DLC step drill $\varnothing 4/\varnothing 4.5$ L10 mm   | 1   | KYFoC6522 |
| Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L11,5 mm | DLC step drill $\varnothing 4/\varnothing 4.5$ L11.5 mm | 1   | KYFoC6523 |

| Descripción - Description                                   |   | Un. | Ref.      |
|---|---|-----|-----------|
| Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L13 mm | DLC step drill $\varnothing 4/\varnothing 4.5$ L13 mm | 1   | KYFoC6524 |
| Fresa cortical DLC $\varnothing 3,2$ mm                     | DLC cortical drill $\varnothing 3.2$ mm               | 1   | KYFoC6401 |
| Fresa cortical DLC $\varnothing 3,55$ mm                    | DLC cortical drill $\varnothing 3.55$ mm              | 1   | KYFoC6402 |
| Fresa cortical DLC $\varnothing 3,9$ mm                     | DLC cortical drill $\varnothing 3.9$ mm               | 1   | KYFoC6403 |
| Fresa cortical DLC $\varnothing 4,4$ mm                     | DLC cortical drill $\varnothing 4.4$ mm               | 1   | KYFoC6404 |
| Fresa cortical DLC $\varnothing 4,9$ mm                     | DLC cortical drill $\varnothing 4.9$ mm               | 1   | KYFoC6405 |
| Llave carraca corta Hex-1,20 mm                             | Short Hex-1.20 mm ratchet wrench                      | 1   | KYLoFo128 |
| Adaptador Hex-4,00 mm LLCIP                                 | Hex-4.00 mm adapter TIRW                              | 1   | KYLoC6054 |
| Adaptador Hex-4,00 mm CA                                    | Hex-4.00 mm adapter HP                                | 1   | KYLoC6055 |
| Extractor transportadores                                   | Mounters extraction wrench                            | 1   | KYLoC6115 |
| Transportadores NP/RP + tornillo                            | Implant carriers NP/RP + screw                        | 2   | KYLoF6107 |
| Transportadores WP + tornillo                               | Implant carriers WP + screw                           | 2   | KYLoF6108 |
| Llave directa NP/RP CA                                      | Direct HP wrench NP/RP                                | 1   | KYLoF6111 |
| Llave directa WP CA   | Direct HP wrench WP                                   | 1   | KYLoF6112 |
| Llave directa NP/RP LLCIP                                   | Direct TIRW NP/RP wrench                              | 1   | KYLoF6113 |
| Llave directa WP LLCIP                                      | Direct TIRW WP wrench                                 | 1   | KYLoF6114 |
| Pin lateral férula  | Lateral pin   | 3   | KYLoC6059 |
| Llave carraca IP  | TI ratchet wrench                                     | 1   | KYLoF0113 |
| Llave transportador plana 45°                               | 45° implant carrier manual wrench                     | 1   | KYLoC0093 |

\*Todas las referencias se venden por separado como recambios  
\*All references are sold separately as spare parts

| Accesorios - Accessories               |                              | Un. | Ref.      |
|--|------------------------------|-----|-----------|
| Casquillo pin lateral férula           | Sleeve for lateral pin       | 5   | KYLoF6061 |
| Casquillo férula $\varnothing 5,10$ mm | Sleeve $\varnothing 5.10$ mm | 5   | KYLoF6102 |

## Kit quirúrgico gmi® PEAK - gmi® PEAK surgical kit

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gmi® PEAK (Ref. KYCoF3054) surgical kit consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of frontier PEAK and avantgard PEAK dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi® surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134°C.

To complement this kit, gmi® has the **frontier/avantgard drill Kit (Ref. KYCoF3055 - see instrumental catalogue)**, allowing it to be easily expanded for the placement of 2 additional implant systems (frontier and avantgard).

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Ref. KYCoF3061 - Kit for gmi® avantgard PEAK system only, without the wrenches exclusive to the frontier PEAK system (see table, highlighted in red).

| Descripción - Description                  |   | Un. | Ref.      |
|--|---|-----|-----------|
| Fresa inicial ø2,0 mm                      | Initial drill ø2.0 mm                       | 1   | KYFoC1401 |
| Fresa piloto ø2,50 mm                      | Pilot drill ø2.50 mm                        | 1   | KYFoC2250 |
| Fresa cónica ø3,20 mm                      | Conical drill ø3.20 mm                      | 1   | KYFoC7032 |
| Fresa cónica ø3,70 mm                      | Conical drill ø3.70 mm                      | 1   | KYFoC7037 |
| Fresa cónica ø4,20 mm                      | Conical drill ø4.20 mm                      | 1   | KYFoC7042 |
| Fresa cortical ø3,20 mm                    | Cortical drill ø3.20 mm                     | 1   | KYFoC4333 |
| Fresa cortical ø3,65 mm                    | Cortical drill ø3.65 mm                     | 1   | KYFoC4337 |
| Fresa cortical ø4,15 mm                    | Cortical drill ø4.15 mm                     | 1   | KYFoC4342 |
| Fresa cortical ø4,65 mm                    | Cortical drill ø4.65 mm                     | 1   | KYFoC4347 |
| Fresa escalonada ø2,50 / ø2,80 mm          | Step drill ø2.50 / ø2.80 mm                 | 1   | KYFoC5129 |
| Fresa escalonada ø3,20 / ø3,65 mm          | Step drill ø3.20 / ø3.65 mm                 | 1   | KYFoC5137 |
| Fresa escalonada ø3,70 / ø4,15 mm          | Step drill ø3.70 / ø4.15 mm                 | 1   | KYFoC5142 |
| Fresa escalonada ø4,15 / ø4,65 mm          | Step drill ø4.15 / ø4.65 mm                 | 1   | KYFoC5147 |
| Indicador profundidad ø2,50 mm             | Depth gauge ø2.50 mm                        | 1   | KYLoC0165 |
| Indicador profundidad cónico ø3,20 mm      | Conical depth gauge ø3.20 mm                | 1   | KYLoC0159 |
| Indicador profundidad cónico ø3,70 mm      | Conical depth gauge ø3.70 mm                | 1   | KYLoC0160 |
| Indicador profundidad cónico ø4,20 mm      | Conical depth gauge ø4.20 mm                | 1   | KYLoC0161 |
| Llave carraca IP                           | TI ratchet wrench                           | 1   | KYLoF0113 |
| Llave-transportador HEX - 2,45 mm CA       | HEX-2.45 mm implant carrier HP wrench       | 1   | KYLoF0091 |
| Llave-transportador HEX - 2,45 mm IP corta | Short HEX-2.45 mm implant carrier TI wrench | 1   | KYLoF0096 |
| Llave-transportador HEX - 2,45 mm IP larga | Long HEX-2.45 mm implant carrier TI wrench  | 1   | KYLoF0111 |
| Llave Avantgard NP/RP CA corta             | Short Avantgard NP/RP HP wrench             | 1   | KYLoF0151 |
| Llave Avantgard NP/RP IP corta             | Short Avantgard NP/RP TI wrench             | 1   | KYLoF0157 |
| Llave Avantgard NP/RP IP larga             | Long Avantgard NP/RP TI wrench              | 1   | KYLoF0162 |
| Alargador fresas                           | Drill extender                              | 1   | KYLoC0023 |
| Paralelizador ø2,00 / ø2,80 mm             | Paralleling pin ø2.00 / ø2.80 mm            | 3   | KYLoC0078 |
| Llave carraca corta HEX-1,20 mm            | Short HEX-1.20 mm ratchet wrench            | 1   | KYLoF0128 |
| Llave carraca larga HEX-1,20 mm            | Long HEX-1.20 mm ratchet wrench             | 1   | KYLoF0129 |

## Kit cirugía guiada gmi® avantgard PEAK - gmi avantgard PEAK guided surgery kit

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gmi® avantgard PEAK (Ref. KYCoF3050) guided surgery kit consists of an autoclavable case made of high strength resin, containing all the necessary material for performing totally guided surgeries with gmi® avantgard PEAK implants of  $\varnothing 3.30 / \varnothing 3.75 / \varnothing 4.25$  and  $\varnothing 4.75$  mm and lengths between 8 and 13 mm.

To complement this kit, gmi® has the frontier/avantgard guided surgery drill Kit (Ref. KYCoF3051 - see instrumental catalogue), allowing it to be easily expanded for the placement of 1 additional implant system (avantgard).

| Descripción - Description                                       |   | Un. | Ref.      |
|---|---|-----|-----------|
| Fresa Pin $\varnothing 1,50$ mm                                 | Pin drill $\varnothing 1.50$ mm                           | 1   | KYFoC6060 |
| Trepante gingival CA $\varnothing 5,10$ mm                      | HP gingival punch $\varnothing 5.10$ mm                   | 1   | KYFoC6066 |
| Fresa inicial $\varnothing 2,5$ mm                              | Initial drill $\varnothing 2.5$ mm                        | 1   | KYFoC6101 |
| Fresa piloto DLC $\varnothing 2,5$ L8 mm                        | DLC pilot drill $\varnothing 2.5$ L8 mm                   | 1   | KYFoC6206 |
| Fresa piloto DLC $\varnothing 2,5$ L10 mm                       | DLC pilot drill $\varnothing 2.5$ L10 mm                  | 1   | KYFoC6207 |
| Fresa piloto DLC $\varnothing 2,5$ L11,5 mm                     | DLC pilot drill $\varnothing 2.5$ L11.5 mm                | 1   | KYFoC6208 |
| Fresa piloto DLC $\varnothing 2,5$ L13 mm                       | DLC pilot drill $\varnothing 2.5$ L13 mm                  | 1   | KYFoC6209 |
| Fresa piloto DLC $\varnothing 2,5$ L15 mm                       | DLC pilot drill $\varnothing 2.5$ L15 mm                  | 1   | KYFoC6210 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L8 mm    | DLC step drill $\varnothing 2.5/\varnothing 2.8$ L8 mm    | 1   | KYFoC6525 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L10 mm   | DLC step drill $\varnothing 2.5/\varnothing 2.8$ L10 mm   | 1   | KYFoC6526 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L11,5 mm | DLC step drill $\varnothing 2.5/\varnothing 2.8$ L11.5 mm | 1   | KYFoC6527 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L13 mm   | DLC step drill $\varnothing 2.5/\varnothing 2.8$ L13 mm   | 1   | KYFoC6528 |
| Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L15 mm   | DLC step drill $\varnothing 2.5/\varnothing 2.8$ L15 mm   | 1   | KYFoC6529 |
| Fresa cónica DLC $\varnothing 3,2$ L8 mm                        | DLC conical drill $\varnothing 3.2$ L8 mm                 | 1   | KYFoC6706 |
| Fresa cónica DLC $\varnothing 3,2$ L10 mm                       | DLC conical drill $\varnothing 3.2$ L10 mm                | 1   | KYFoC6707 |
| Fresa cónica DLC $\varnothing 3,2$ L11,5 mm                     | DLC conical drill $\varnothing 3.2$ L11.5 mm              | 1   | KYFoC6708 |
| Fresa cónica DLC $\varnothing 3,2$ L13 mm                       | DLC conical drill $\varnothing 3.2$ L13 mm                | 1   | KYFoC6709 |
| Fresa cónica DLC $\varnothing 3,2$ L15 mm                       | DLC conical drill $\varnothing 3.2$ L15 mm                | 1   | KYFoC6710 |
| Fresa cónica DLC $\varnothing 3,7$ L8 mm                        | DLC conical drill $\varnothing 3.7$ L8 mm                 | 1   | KYFoC6711 |
| Fresa cónica DLC $\varnothing 3,7$ L10 mm                       | DLC conical drill $\varnothing 3.7$ L10 mm                | 1   | KYFoC6712 |
| Fresa cónica DLC $\varnothing 3,7$ L11,5 mm                     | DLC conical drill $\varnothing 3.7$ L11.5 mm              | 1   | KYFoC6713 |
| Fresa cónica DLC $\varnothing 3,7$ L13 mm                       | DLC conical drill $\varnothing 3.7$ L13 mm                | 1   | KYFoC6714 |
| Fresa cónica DLC $\varnothing 3,7$ L15 mm                       | DLC conical drill $\varnothing 3.7$ L15 mm                | 1   | KYFoC6715 |
| Fresa cónica DLC $\varnothing 4,2$ L8 mm                        | DLC conical drill $\varnothing 4.2$ L8 mm                 | 1   | KYFoC6716 |
| Fresa cónica DLC $\varnothing 4,2$ L10 mm                       | DLC conical drill $\varnothing 4.2$ L10 mm                | 1   | KYFoC6717 |
| Fresa cónica DLC $\varnothing 4,2$ L11,5 mm                     | DLC conical drill $\varnothing 4.2$ L11.5 mm              | 1   | KYFoC6718 |

| Descripción - Description                 |  | Un. | Ref.      |
|---|--|-----|-----------|
| Fresa cónica DLC $\varnothing 4,2$ L13 mm | DLC conical drill $\varnothing 4.2$ L13 mm | 1   | KYFoC6719 |
| Fresa cortical DLC $\varnothing 3,2$ mm   | DLC cortical drill $\varnothing 3.2$ mm    | 1   | KYFoC6406 |
| Fresa cortical DLC $\varnothing 3,65$ mm  | DLC cortical drill $\varnothing 3.65$ mm   | 1   | KYFoC6407 |
| Fresa cortical DLC $\varnothing 4,15$ mm  | DLC cortical drill $\varnothing 4.15$ mm   | 1   | KYFoC6408 |
| Fresa cortical DLC $\varnothing 4,65$ mm  | DLC cortical drill $\varnothing 4.65$ mm   | 1   | KYFoC6409 |
| Llave carraca corta Hex-1,20 mm           | Short Hex-1.20 mm ratchet wrench           | 1   | KYLoF0128 |
| Adaptador Hex-4,00 mm LLCIP               | Hex-4.00 mm adapter TIRW                   | 1   | KYLoC6054 |
| Adaptador Hex-4,00 mm CA                  | Hex-4.00 mm adapter HP                     | 1   | KYLoC6055 |
| Extractor transportadores                 | Mounters extraction wrench                 | 1   | KYLoC6115 |
| Transportadores NP/RP + tornillo          | Implant carriers NP/RP + screw             | 2   | KYLoF6107 |
| Transportadores WP + tornillo             | Implant carriers WP + screw                | 2   | KYLoF6108 |
| Llave directa NP/RP CA                    | Direct HP wrench NP/RP                     | 1   | KYLoF6111 |
| Llave directa WP CA                       | Direct HP wrench WP                        | 1   | KYLoF6112 |
| Llave directa NP/RP LLCIP                 | Direct TIRW NP/RP wrench                   | 1   | KYLoF6113 |
| Llave directa WP LLCIP                    | Direct TIRW WP wrench                      | 1   | KYLoF6114 |
| Pin lateral férula                        | Lateral pin                                | 3   | KYLoC6059 |
| Llave carraca IP                          | TI ratchet wrench                          | 1   | KYLoF0113 |
| Llave transportador plana 45°             | 45° implant carrier manual wrench          | 1   | KYLoC0093 |

\*Todas las referencias se venden por separado como recambios  
 \*All references are sold separately as spare parts

| Accesorios - Accessories               |                              |     |           |
|--|------------------------------|-----|-----------|
| Descripción - Description              |                              | Un. | Ref.      |
| Casquillo pin lateral férula           | Sleeve for lateral pin       | 5   | KYLoF6061 |
| Casquillo férula $\varnothing 5,10$ mm | Sleeve $\varnothing 5.10$ mm | 5   | KYLoF6102 |

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