# SISTEMAS DE IMPLANTES IMPLANT SYSTEMS

www.gmidental.com

avantgard

avantgard  
PEAK

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### KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

### IMPLANT SYSTEMS

- **gmi avantgard** dental implants 8
- Insertion sequence for **gmi avantgard** implants 9
- **gmi avantgard** implants range 10
- Drilling sequences for **gmi avantgard** implants 12
- **gmi avantgard PEAK** dental implants 14
- Insertion sequence for **gmi avantgard PEAK** implants 15
- **gmi avantgard PEAK** implants range 16
- Drilling sequences for **gmi avantgard PEAK** implants 17

### ATTACHMENTS

- **gmi avantgard** healing abutments 19
- **gmi avantgard** implant carrier abutments 20
- **gmi avantgard** temporary cylinders 21
- **gmi avantgard** angled abutments 22
- **gmi avantgard** straight abutments 23
- **gmi avantgard** titanium base 24
- **gmi avantgard** multi-esthetic abutments 26
- **gmi avantgard** bioesthetic bases 29
- **gmi avantgard** CrCo pre-milled blanks 30
- **gmi avantgard** replicas 31
- **gmi avantgard** impression coping 32
- **gmi avantgard** clinic and laboratory screws 33
- **gmi avantgard** MEA and BE screws 34
- Recommended tightening torques 35

### GMI SURGICAL KITS

- **gmi frontier/avantgard** surgical kit 37
- **gmi avantgard** guided surgery kit 38
- **gmi PEAK** surgical kit 39
- **gmi avantgard PEAK** guided surgery kit 40

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


To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.





SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

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The **gmi avantgard** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

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The **gmi avantgard** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi avantgard** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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


To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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
In the chart are detailed all the symbols appearing on the gmi avantgard implant system's labelling and packaging and their corresponding description.

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\* Minimum requirements for downloading and opening the electronic instructions for use:  
- Stable internet connectivity.  
- A device with a modern operating system that supports PDF viewing.  
- Adequate storage space on your device.  
- A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.  
- A PDF reader application installed on your device..

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
gmi avantgard conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high primary stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

### Características generales Key features

- |   |  |
|---|--|
| <p>① Plataformas<br/>NP: estrecha<br/>RP: estándar<br/>WP: ancha</p> <p>② Longitud</p> <p>③ Fresados auto-roscales<br/>Longitud variable<br/>3 fresados a 120°<br/>Ápice atraumático</p> <p>④ Rosca cortical</p> <p>⑤ Regatas internas<br/>6 x 60°</p> <p>⑥ Rosca interior<br/>NP-RP: M1,6 x 0,35 mm<br/>WP: M2,0 x 0,40mm</p> <p>⑦ Cuerpo del implante<br/>Paso: 1,60 mm</p> | <p>① Platforms<br/>NP: narrow<br/>RP: standard<br/>WP: wide</p> <p>② Length</p> <p>③ Self threading millings<br/>Variable length<br/>3 millings at 120°<br/>Atraumatic tip</p> <p>④ Cortical thread</p> <p>⑤ Internal grooves<br/>6 x 60°</p> <p>⑥ Internal thread<br/>NP-RP: M1.6 x 1.35 mm<br/>WP: M2.0 x 0.40 mm</p> <p>⑦ Implant body<br/>Pitch: 1.60 mm</p> |
|---|--|

2797

	NP	RP	RP	WP	WP
①	NP	RP	RP	WP	WP
②	De 10 a 15 From 10 to 15	De 8 a 15 From 8 to 15	De 6,5 a 13 From 6.5 to 13	De 6,5 a 13 From 6.5 to 13	De 6,5 a 11,5 From 6.5 to 11.5
④	ø3,30	ø4,00	ø4,50	ø5,00	ø6,00
⑦	ø3,30	ø3,75	ø4,25	ø4,75	ø5,75

## Protocolo de inserción de implantes gmi avantgard - Insertion sequence for gmi avantgard implants

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**gmi avantgard** conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi avantgard** implants consists of a titanium support which houses the dental implant, prepared to be inserted directly into the mouth, using the implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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\* Consult the instructions for use in electronic format.


1. Check sticker indicating implant diameter and length and open the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the insertion wrench facing the grooves of the wrench with those of the implant.
7. Check that the wrench is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the subcrestal level. Remove the wrench by making an axial movement.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove wrench and suture the incision.

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


The following is a description of the drilling sequences for the different conical connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped drill: 1200 - 1500 rpm.
- Pilot drill: 700 - 900 rpm.
- Step drills:
  - Ø3.00 & Ø3.50 mm → 400 - 700 rpm.
  - Ø4.00 & Ø5.10 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.

For drilling, external cooling with saline solution is required.


avantgard Ø3,30			
Descripción - Description			Ref.
1	Fresa lanza	Lance-shaped drill	KYFoC1225
2	Fresa piloto Ø2,50	Pilot drill Ø2.50	KYFoC2250
3	Fresa cortical Ø3,20	Cortical drill Ø3,20	KYFoC4332
4	Fresa escalonada Ø2,50/Ø3,00	Step drill Ø2.50/Ø3.00	KYFoC5130


avantgard Ø3,75			
Descripción - Description			Ref.
1	Fresa lanza	Lance-shaped drill	KYFoC1225
2	Fresa piloto Ø2,50	Pilot drill Ø2.50	KYFoC2250
3	Fresa escalonada Ø2,50/Ø3,00	Step drill Ø2.50/Ø3.00	KYFoC5130
4	Fresa cortical Ø4,00	Cortical drill Ø4.00	KYFoC4340
5	Fresa escalonada Ø3,00/Ø3,50	Step drill Ø3.00/Ø3.50	KYFoC5135

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

gmi avantgard PEAK conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the conical core and progressive thread profile facilitate implant insertion and act as a bone compactor allowing high primary stability to be obtained in sparse bones.

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**gmi avantgard PEAK** conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi avantgard PEAK** implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single Avantgard NP-RP implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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\* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the carton.
2. **Remove the secondary packaging and product identification labels.**
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. **Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.**
5. Extract and save the primary packaging plug making a lateral movement.
6. **Insert the NP-RP insertion wrench facing the grooves of the wrench with those of the implant.**
7. Check that the wrench is fully inserted and turn slightly while gently pulling on the assembly.
8. **Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the subcrestal level. Remove the wrench by making an axial movement.**
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. **Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove wrench and suture the incision.**

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

The following is a description of the drilling sequences for the different conical connection implant models, as well as the recommended drilling conditions for their use:

- Initial drill: 1000 rpm.
- Pilot drill: 700 - 900 rpm.
- Step and conical drills:
  - ∅2.50 / ∅2.80 mm → 500 - 700 rpm.
  - ∅3.20 / ∅3.65 mm → 500 - 700 rpm.
  - ∅3.70 / ∅4.15 mm → 400 - 700 rpm.
  - ∅4.15 / ∅4.65 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.
- For drilling, external cooling with saline solution is required.


avantgard PEAK Ø3,30			
Descripción - Description		Ref.	
1	Fresa inicial	Initial drill	KYFOC1401
2	Fresa piloto Ø2,50 mm	Pilot drill Ø2.50	KYFOC2250
3	Fresa escalonada Ø2,50/Ø2,80 mm	Step drill Ø2.50/Ø2.80 mm	KYFOC5129
4	Fresa cortical Ø3,20 mm	Cortical drill Ø3.20	KYFOC4333


avantgard PEAK Ø3,75			
Descripción - Description		Ref.	
1	Fresa inicial	Initial drill	KYFOC1401
2	Fresa piloto Ø2,50 mm	Pilot drill Ø2.50	KYFOC2250
3	Fresa escalonada Ø2,50/Ø2,80 mm	Step drill Ø2.00/Ø2.80 mm	KYFOC5129
4	Fresa cónica Ø3,20 mm	Conical drill Ø3.20 mm	KYFOC7032
5	Fresa cortical Ø3,65 mm	Cortical drill Ø3.65	KYFOC4337
6	Fresa escalonada Ø3,20/Ø3,65 mm	Step drill Ø3.20/Ø3.65 mm	KYFOC5137

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Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Available in emergency heights from 1.50 to 4.50 mm.
- Apical bit to facilitate the start of threading.
- Manual tightening torque: Max. 15 N-cm.

Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
NP ●●	Ø 4,00	1,50	KDPoC4012
		2,50	KDPoC4013
		3,50	KDPoC4014
		4,50	KDPoC4015
RP ●●	Ø 5,00	1,50	KDPoC4016
		2,50	KDPoC4017
		3,50	KDPoC4018
		4,50	KDPoC4019

Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
NP ●●	Ø 5,00	1,50	KDPoC4001
		2,50	KDPoC4002
		3,50	KDPoC4003
		4,50	KDPoC4004
RP ●●	Ø 6,00	2,50	KDPoC4009
		3,50	KDPoC4010
		4,50	KDPoC4011

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
When the use of a custom abutment is required, implant carrier abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, implant carrier abutments can be used like closed tray impression coping.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 3.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping.
- Premounted clinic screw included.
- Screw tightening torque: 25 N-cm.

Plataforma Platform	Ø hombro Ø shoulder	Ref.
RP ●●	Ø 3,75 - Ø 4,25	KDBoF4001
WP ●	Ø 4,75	KDBoF4002
	Ø 5,75	KDBoF4003

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
When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary abutment is used, defined as attachment that once attached to the crown of the implant with clinic screw, they will allow a temporary restoration to be carried out while the definitive restoration is being manufactured.

### Key features

- Machined in high strength technical plastic.
- NP abutments machined in Grade V titanium as per ISO-5832 standard
- Rotary version only for multiple restorations.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw included.
- Manual tightening torque: Max. 15 N·cm.

Plataforma Platform	Tipo conexión Connection type	Ref.
RP ● ●	Anti-rotatoria • Non-rotary	KDHOF4001
	Rotatoria • Rotary	KDHOF4002
WP ●	Anti-rotatoria • Non-rotary	KDHOF4003
	Rotatoria • Rotary	KDHOF4004

Plataforma Platform	Altura (mm) Height (mm)	Tipo conexión Connection type	Ref.
NP ● ●	1,50	Anti-rotatoria • Non-rotary	KDHOF4009
	3,00	Anti-rotatoria • Non-rotary	KDHOF4010
	1,50	Rotatoria • Rotary	KDHOF4011
	3,00	Rotatoria • Rotary	KDHOF4012

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw included.
- Abutment shoulder adapted to the shape of the gum.
- Screw tightening torque: 25 N-cm.

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
NP ●●	10°	2,50	KDMoF4039
	10°	4,50	KDMoF4040
	20°	2,50	KDMoF4041
	20°	4,50	KDMoF4042

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP ●●	15°	2,50	KDMoF4013
	15°	4,50	KDMoF4014
	20°	2,50	KDMoF4015
	20°	4,50	KDMoF4016

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
WP ●	15°	2,50	KDMoF4017
	15°	4,50	KDMoF4018
	20°	2,50	KDMoF4019
	20°	4,50	KDMoF4020

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
Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

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
### Key features


- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Available only in non-rotary version.
- Clinic screw included.
- Available in four emergency heights from 1.50 to 4.50 mm.
- Screw tightening torque: 25 N-cm.

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
NP ●●	Ø 4,00	1,50	KDMoF4035
		2,50	KDMoF4036
		3,50	KDMoF4037
		4,50	KDMoF4038

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
RP ●●	Ø 4,00	1,50	KDMoF4001
		2,50	KDMoF4002
		3,50	KDMoF4003
		4,50	KDMoF4004
	Ø 5,00	1,50	KDMoF4005
		2,50	KDMoF4006
		3,50	KDMoF4007
		4,50	KDMoF4008

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
WP ●	Ø 6,00	1,50	KDMoF4009
		2,50	KDMoF4010
		3,50	KDMoF4011
		4,50	KDMoF4012

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Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw included.
- Available in emergency height of 1.50 mm and 3.00 mm.
- Non-rotary titanium base with indexation.
- Screw tightening torque: 25 N-cm.
- Scan directly on Titanium base or alternative scan on Scan abutment.

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

Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.50 to 4.50 mm.
- One-piece straight abutments in rotary version.
- Tightened with hex wrench: B/F 2.00 mm (KYLOCO149).
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.
- Tightening torques: ME straight abut. (30 N-cm), ME angled abut. (20 N-cm) and healing caps (Manual: Max. 15 N-cm).

### Pilares multi-estético rectos Multi-esthetic straight abutments

Plataforma Platform	Altura (mm) Height (mm)	Ref.
NP ●●	1,50	KDGoF4019
	2,50	KDGoF4020
	3,50	KDGoF4021
	4,50	KDGoF4022
RP ●●	1,50	KDGoF4001
	2,50	KDGoF4002
	3,50	KDGoF4003
	4,50	KDGoF4004
WP ●	1,50	KDGoF4005
	2,50	KDGoF4006
	3,50	KDGoF4007
	4,50	KDGoF4008

## Pilares multi-estético gmi avantgard - gmi avantgard multi-esthetic abutments

### Pilares multi-estético inclinados Multi-esthetic angled abutments

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
	17°	2,50	KDGoF4009
RP ●●	17°	3,50	KDGoF4010
	30°	4,00	KDGoF4011


### Pilares cicatrización PME Healing abutments MEA

Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
NP - RP - WP ●●	Ø 4,80	KDPoF0035
	Ø 6,00	KDPoF0037

### Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

Plataforma Platform	Ref.
NP - RP - WP ●●	KDSoF0015

### Tomas impresión cubeta cerrada (CC) PME Closed tray (CT) impression coping MEA

Plataforma Platform	Cantidad (Ud.) Quantity (Un.)	Ref.
	1	KDSoF0017
NP - RP - WP ●●	5	KDSoF5001

## Pilares multi-estético gmi avantgard - gmi avantgard multi-esthetic abutments

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Once the insertion of the implant in the bone bed has been completed, the bioesthetic base is used, which is defined as the machined attachment that, once fixed to the coronal part of the implant, allows single restorations at the abutment level.

Placing the bioesthetic base during surgery allows the soft tissues to heal around it.

Once the tissue has healed, it serves as the basis for the rehabilitation, avoiding the connection and disconnection of the components to the implant and, therefore, keeping the biological seal intact.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 2.50 to 4.50 mm.
- Straight two-piece bases, including the screw.
- Tightened with hex wrench: B/F 2.00 mm (KYLOCo149).
- Impression coping with screw included.
- Castable.
- Tightening torques: BE base (25 N-cm) and healing abutment (Manual: Max. 15 N-cm).

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
CrCo pre-milled blanks consist of a block with the implant connection machined, which allows an excellent fit between the implant and the abutment.

They are indicated to produce custom one-piece abutments using CAD/CAM allowing direct application of ceramic.

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

Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.

Réplicas de implante Implant replicas	
Plataforma Platform	Ref.
NP ●●	KDRoC4005
RP ●●	KDRoC4001
WP ●	KDRoC4002

Réplicas digitales Digital replicas	
Plataforma Platform	Ref.
NP ●●	KDRoF4006
RP ●●	KDRoF4003
WP ●	KDRoF4004

Réplicas para pilares multi-estético MEA replicas	
Plataforma Platform	Ref.
NP - RP - WP ●●	KDRoC0012

Réplicas digitales para pilares multi-estético MEA digital replicas	
Plataforma Platform	Ref.
NP - RP - WP ●●	KDRoF0014

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
Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, is used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

### Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- Manual tightening torque: Max. 15 N-cm.

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
Plataforma Platform	Ref.
NP ●●	KDSOF4007
	KDSOF0019 (x5)
RP - WP ●●	KDSOF0014 (x5)

## Tornillos clínica y laboratorio gmi avantgard - gmi avantgard clinic and laboratory screws

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
Clinic screws are defined as machined devices used for the final anchoring of the abutments or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of abutments to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

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### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo clínica Clinic screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
NP - RP ●●	⊕ 1,20 mm	KDToC4003
	⊕ TX	KDToC4021
WP ●	⊕ 1,20 mm	KDToC4004
	⊕ TX	KDToC4022
Tornillo laboratorio Laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
NP - RP ●●	⊕ 1,20 mm	KDToC4005
	⊕ TX	KDToC4023
WP ●	⊕ 1,20 mm	KDToC4006
	⊕ TX	KDToC4024

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MEA and BE clinical screws are defined as machined devices used for the final anchoring of the ME and BE abutments, while laboratory screws are the machined devices used for the temporary anchoring of ME and BE abutments to the replica during the fabrication of the prostheses in the laboratory.

### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm and 2.00 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo PME MEA screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
NP - RP - WP ● ●	○ 1,20 mm	KDToC0056
	⊛ TX	KDToC0070

Tornillo laboratorio PME MEA laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
NP - RP - WP ● ●	○ 1,20 mm	KDToC0057
	⊛ TX	KDToC0072

Tornillo BE BE screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ● ●	○ 1,20 mm	KDToC4016
	⊛ TX	KDToC0074

Tornillo laboratorio BE BE laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ● ●	○ 1,20 mm	KDToC4020
	⊛ TX	KDToC0075

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

**gmi frontier/avantgard surgical kit** consist of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi frontier and avantgard dental implant systems**. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

**gmi surgical kits** are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134 °C.

To complement this kit, **gmi** has the **PEAK drill Kit** (*see instrumental catalogue*), allowing it to be easily expanded for the placement of 2 additional implant systems (frontier PEAK and avantgard PEAK).

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gmi avantgard guided surgery kit consists of an autoclavable case made of high strength resin, containing all the necessary material for performing totally guided surgeries with gmi avantgard implants of  $\varnothing 3.30$  /  $\varnothing 3.75$  /  $\varnothing 4.25$  and  $\varnothing 4.75$  mm and lengths between 8 and 13 mm.

To complement this kit, gmi has the PEAK guided surgery drill Kit (see instrumental catalogue), allowing it to be easily expanded for the placement of 1 additional implant system (avantgard PEAK).

Descripción - Description		Un.	Ref.
Fresa Pin $\varnothing 1,50$ mm	Pin drill $\varnothing 1.50$ mm	1	KYFOC6060
Trepante gingival CA $\varnothing 5,10$ mm	HP gingival punch $\varnothing 5.10$ mm	1	KYFOC6066
Fresa inicial $\varnothing 2,5$ mm	Initial drill $\varnothing 2.5$ mm	1	KYFOC6101
Fresa piloto DLC $\varnothing 2,5$ L8 mm	DLC pilot drill $\varnothing 2.5$ L8 mm	1	KYFOC6201
Fresa piloto DLC $\varnothing 2,5$ L10 mm	DLC pilot drill $\varnothing 2.5$ L10 mm	1	KYFOC6202
Fresa piloto DLC $\varnothing 2,5$ L11,5 mm	DLC pilot drill $\varnothing 2.5$ L11.5 mm	1	KYFOC6203
Fresa piloto DLC $\varnothing 2,5$ L13, mm	DLC pilot drill $\varnothing 2.5$ L13 mm	1	KYFOC6204
Fresa piloto DLC $\varnothing 2,5$ L15 mm	DLC pilot drill $\varnothing 2.5$ L15 mm	1	KYFOC6205
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L8 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L8 mm	1	KYFOC6506
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L10 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L10 mm	1	KYFOC6507
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L11,5 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L11.5 mm	1	KYFOC6508
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L13 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L13 mm	1	KYFOC6509
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L15 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L15 mm	1	KYFOC6510
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L8 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L8 mm	1	KYFOC6511
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L10 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L10 mm	1	KYFOC6512
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L11,5 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L11.5 mm	1	KYFOC6513
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L13 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L13 mm	1	KYFOC6514
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L15 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L15 mm	1	KYFOC6515
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L8 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L8 mm	1	KYFOC6516
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L10 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L10 mm	1	KYFOC6517
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L11,5 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L11.5 mm	1	KYFOC6518
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L13 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L13 mm	1	KYFOC6519
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L15 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L15 mm	1	KYFOC6520
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L8 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L8 mm	1	KYFOC6521
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L10 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L10 mm	1	KYFOC6522
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L11,5 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L11.5 mm	1	KYFOC6523

Descripción - Description		Un.	Ref.
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L13 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L13 mm	1	KYFOC6524
Fresa cortical DLC $\varnothing 3,2$ mm	DLC cortical drill $\varnothing 3.2$ mm	1	KYFOC6401
Fresa cortical DLC $\varnothing 3,55$ mm	DLC cortical drill $\varnothing 3.55$ mm	1	KYFOC6402
Fresa cortical DLC $\varnothing 3,9$ mm	DLC cortical drill $\varnothing 3.9$ mm	1	KYFOC6403
Fresa cortical DLC $\varnothing 4,4$ mm	DLC cortical drill $\varnothing 4.4$ mm	1	KYFOC6404
Fresa cortical DLC $\varnothing 4,9$ mm	DLC cortical drill $\varnothing 4.9$ mm	1	KYFOC6405
Llave carraca corta Hex-1,20 mm	Short Hex-1.20 mm ratchet wrench	1	KYLOF0128
Adaptador Hex-4,00 mm LLCIP	Hex-4.00 mm adapter TIRW	1	KYLOC6054
Adaptador Hex-4,00 mm CA	Hex-4.00 mm adapter HP	1	KYLOC6055
Extractor transportadores	Mounters extraction wrench	1	KYLOC6115
Transportadores NP/RP + tornillo	Implant carriers NP/RP + screw	2	KYLOF6107
Transportadores WP + tornillo	Implant carriers WP + screw	2	KYLOF6108
Llave directa NP/RP CA	Direct HP wrench NP/RP	1	KYLOF6111
Llave directa WP CA	Direct HP wrench WP	1	KYLOF6112
Llave directa NP/RP LLCIP	Direct TIRW NP/RP wrench	1	KYLOF6113
Llave directa WP LLCIP	Direct TIRW WP wrench	1	KYLOF6114
Pin lateral férula	Lateral pin	3	KYLOF6059
Llave carraca IP	TI ratchet wrench	1	KYLOF0113
Llave transportador plana 45°	45° implant carrier manual wrench	1	KYLOC0093

\*Todas las referencias se venden por separado como recambios  
\*All references are sold separately as spare parts

Accesorios - Accessories		Un.	Ref.
Casquillo pin lateral férula	Sleeve for lateral pin	5	KYLOF6061
Casquillo férula $\varnothing 5,10$ mm	Sleeve $\varnothing 5.10$ mm	5	KYLOF6102

## Kit quirúrgico gmi PEAK - gmi PEAK surgical kit

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**gmi PEAK surgical kit** consist of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi frontier PEAK** and **avantgard PEAK** dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

**gmi surgical kits** are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134 °C.

To complement this kit, **gmi** has the **frontier/avantgard drill Kit** (see *instrumental catalogue*), allowing it to be easily expanded for the placement of 2 additional implant systems (frontier and avantgard).

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gmi avantgard PEAK guided surgery kit consists of an autoclavable case made of high strength resin, containing all the necessary material for performing totally guided surgeries with gmi avantgard PEAK implants of  $\varnothing 3.30$  /  $\varnothing 3.75$  /  $\varnothing 4.25$  and  $\varnothing 4.75$  mm and lengths between 8 and 13 mm.

To complement this kit, gmi has the **frontier/avantgard guided surgery drill Kit** (see instrumental catalogue), allowing it to be easily expanded for the placement of 1 additional implant system (avantgard).

Descripción - Description		Un.	Ref.
Fresa Pin $\varnothing 1,50$ mm	Pin drill $\varnothing 1.50$ mm	1	KYFoC6060
Trepante gingival CA $\varnothing 5,10$ mm	HP gingival punch $\varnothing 5.10$ mm	1	KYFoC6066
Fresa inicial $\varnothing 2,5$ mm	Initial drill $\varnothing 2.5$ mm	1	KYFoC6101
Fresa piloto DLC $\varnothing 2,5$ L8 mm	DLC pilot drill $\varnothing 2.5$ L8 mm	1	KYFoC6206
Fresa piloto DLC $\varnothing 2,5$ L10 mm	DLC pilot drill $\varnothing 2.5$ L10 mm	1	KYFoC6207
Fresa piloto DLC $\varnothing 2,5$ L11,5 mm	DLC pilot drill $\varnothing 2.5$ L11.5 mm	1	KYFoC6208
Fresa piloto DLC $\varnothing 2,5$ L13 mm	DLC pilot drill $\varnothing 2.5$ L13 mm	1	KYFoC6209
Fresa piloto DLC $\varnothing 2,5$ L15 mm	DLC pilot drill $\varnothing 2.5$ L15 mm	1	KYFoC6210
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L8 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L8 mm	1	KYFoC6525
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L10 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L10 mm	1	KYFoC6526
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L11,5 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L11.5 mm	1	KYFoC6527
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L13 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L13 mm	1	KYFoC6528
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L15 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L15 mm	1	KYFoC6529
Fresa cónica DLC $\varnothing 3,2$ L8 mm	DLC conical drill $\varnothing 3.2$ L8 mm	1	KYFoC6706
Fresa cónica DLC $\varnothing 3,2$ L10 mm	DLC conical drill $\varnothing 3.2$ L10 mm	1	KYFoC6707
Fresa cónica DLC $\varnothing 3,2$ L11,5 mm	DLC conical drill $\varnothing 3.2$ L11.5 mm	1	KYFoC6708
Fresa cónica DLC $\varnothing 3,2$ L13 mm	DLC conical drill $\varnothing 3.2$ L13 mm	1	KYFoC6709
Fresa cónica DLC $\varnothing 3,2$ L15 mm	DLC conical drill $\varnothing 3.2$ L15 mm	1	KYFoC6710
Fresa cónica DLC $\varnothing 3,7$ L8 mm	DLC conical drill $\varnothing 3.7$ L8 mm	1	KYFoC6711
Fresa cónica DLC $\varnothing 3,7$ L10 mm	DLC conical drill $\varnothing 3.7$ L10 mm	1	KYFoC6712
Fresa cónica DLC $\varnothing 3,7$ L11,5 mm	DLC conical drill $\varnothing 3.7$ L11.5 mm	1	KYFoC6713
Fresa cónica DLC $\varnothing 3,7$ L13 mm	DLC conical drill $\varnothing 3.7$ L13 mm	1	KYFoC6714
Fresa cónica DLC $\varnothing 3,7$ L15 mm	DLC conical drill $\varnothing 3.7$ L15 mm	1	KYFoC6715
Fresa cónica DLC $\varnothing 4,2$ L8 mm	DLC conical drill $\varnothing 4.2$ L8 mm	1	KYFoC6716
Fresa cónica DLC $\varnothing 4,2$ L10 mm	DLC conical drill $\varnothing 4.2$ L10 mm	1	KYFoC6717
Fresa cónica DLC $\varnothing 4,2$ L11,5 mm	DLC conical drill $\varnothing 4.2$ L11.5 mm	1	KYFoC6718

Descripción - Description		Un.	Ref.
Fresa cónica DLC $\varnothing 4,2$ L13 mm	DLC conical drill $\varnothing 4.2$ L13 mm	1	KYFoC6719
Fresa cortical DLC $\varnothing 3,2$ mm	DLC cortical drill $\varnothing 3.2$ mm	1	KYFoC6406
Fresa cortical DLC $\varnothing 3,65$ mm	DLC cortical drill $\varnothing 3.65$ mm	1	KYFoC6407
Fresa cortical DLC $\varnothing 4,15$ mm	DLC cortical drill $\varnothing 4.15$ mm	1	KYFoC6408
Fresa cortical DLC $\varnothing 4,65$ mm	DLC cortical drill $\varnothing 4.65$ mm	1	KYFoC6409
Llave carraca corta Hex-1,20 mm	Short Hex-1.20 mm ratchet wrench	1	KYLoF0128
Adaptador Hex-4,00 mm LLCIP	Hex-4.00 mm adapter TIRW	1	KYLoC6054
Adaptador Hex-4,00 mm CA	Hex-4.00 mm adapter HP	1	KYLoC6055
Extractor transportadores	Mounters extraction wrench	1	KYLoC6115
Transportadores NP/RP + tornillo	Implant carriers NP/RP + screw	2	KYLoF6107
Transportadores WP + tornillo	Implant carriers WP + screw	2	KYLoF6108
Llave directa NP/RP CA	Direct HP wrench NP/RP	1	KYLoF6111
Llave directa WP CA	Direct HP wrench WP	1	KYLoF6112
Llave directa NP/RP LLCIP	Direct TIRW NP/RP wrench	1	KYLoF6113
Llave directa WP LLCIP	Direct TIRW WP wrench	1	KYLoF6114
Pin lateral férula	Lateral pin	3	KYLoC6059
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave transportador plana 45°	45° implant carrier manual wrench	1	KYLoC0093

\*Todas las referencias se venden por separado como recambios  
 \*All references are sold separately as spare parts

Accesorios - Accessories			
Descripción - Description		Un.	Ref.
Casquillo pin lateral férula	Sleeve for lateral pin	5	KYLoF6061
Casquillo férula $\varnothing 5,10$ mm	Sleeve $\varnothing 5.10$ mm	5	KYLoF6102

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