## SISTEMAS DE IMPLANTES IMPLANT SYSTEMS

[www.gmidental.com](http://www.gmidental.com)

Implantes dentales de conexión interna hexagonal  
Hexagonal internal connection dental implants

Made in Spain

## Índice - Summary

### CARACTERÍSTICAS GENERALES

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### KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

### IMPLANT SYSTEMS

- **gmi® frontier** dental implants 8
- Insertion sequence for **gmi® frontier** implants 9
- **gmi® frontier** implants range 10
- Drilling sequences for **gmi® frontier** implants 12
- **gmi® frontier PEAK** dental implants 14
- Insertion sequence for **gmi® frontier PEAK** implants 15
- **gmi® frontier PEAK** implants range 16
- Drilling sequences for **gmi® frontier PEAK** implants 17

### GMI® FRONTIER ATTACHMENTS

- Healing abutments 18
- Millable-implant abutments 19
- Esthetic angled abutments 20
- Straight abutments 20
- Temporary cylinders 21
- Titanium base 22
- Ball abutments 23
- Multi-esthetic abutments 24
- Bioesthetic bases 27
- Replicas 29
- Impression coping 30
- Castables 31
- CrCo bases 32
- Clinic and laboratory screws 32
- MEA and BE screws 33
- Recommended tightening torques 35

### GMI® SURGICAL KITS

- **gmi® frontier/avantgard** surgical kit 37
- **gmi® frontier** guided surgery kit 38
- **gmi® PEAK** surgical kit 39
- **gmi® frontier PEAK** guided surgery kit 40

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- The **CE 2797** marking applies to medical devices certified through conformity assessment by a Notified Body under Regulation (EU) 2017/745.

- The **CE** marking (without an identification number) applies to self-certified medical devices under Regulation (EU) 2017/745.




- Neither of these symbols applies to products that fall outside the scope of Regulation (EU) 2017/745 for medical devices.

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To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.








SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

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The **gmi**® **frontier** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

		
<b>ø3,30-ø3,75 mm</b>	<b>ø4,25 mm</b>	<b>ø4,75-ø5,75 mm</b>
<b>Incisivos</b> (Excepto centrales maxilares) Incisors (Except central maxillary incisors)	<b>Caninos · Canines</b> <b>Premolares · Premolars</b> <b>Incisivos centrales maxilares</b> Central maxillary incisors	<b>Molares</b> Molars

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The **gmi**® **frontier** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi**® recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi**® **frontier** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi**® recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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
To prevent prior deformation of the clinic screw, **gmi**® recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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
In the chart are detailed all the symbols appearing on the **gmi® frontier** implant system's labelling and packaging and their corresponding description.

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\* **Minimum requirements for downloading and opening the electronic instructions for use:**  
 - Stable internet connectivity.  
 - A device with a modern operating system that supports PDF viewing.  
 - Adequate storage space on your device.  
 - A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.  
 - A PDF reader application installed on your device.

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
gmi® frontier internal connection hexagonal implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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- ① Platforms  
RP: Standard ø3.30 mm  
WP: Wide ø4.30 mm
- ② Length
- ③ Self threading millings  
Variable length  
3 millings at 120°  
Atraumatic tip
- ④ Cortical thread
- ⑤ Internal hexagon  
B/F 2.45 mm x 1.50 mm
- ⑥ Internal thread  
M1.8 x 0.35 mm
- ⑦ Implant body  
Pitch: 1.60 mm

	●	●	●	●	●
①	RP	RP	RP	WP	WP
②	De 8 a 15 From 8 to 15	De 8 a 17 From 8 to 17	De 6,5 a 17 From 6.5 to 17	De 6,5 a 13 From 6.5 to 13	De 6,5 a 13 From 6.5 to 13
④	ø3,60	ø4,00	ø4,50	ø5,00	ø6,00
⑦	ø3,30	ø3,75	ø4,25	ø4,75	ø5,75

## Protocolo de inserción de implantes gmi® frontier - Insertion sequence for gmi® frontier implants

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gmi® frontier internal connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the gmi® frontier implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single hex-2.45 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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\* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the HEX-2.45 mm implant carrier wrench until a slight retention is noticed, addressing the hex of the key with the faces of the implant hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the crestal or subcrestal level. Remove the key.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove key and suture the incision.

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

The following is a description of the drilling sequences for the different internal connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped drill: 1200 - 1500 rpm.
- Pilot drill: 700 - 900 rpm.
- Step drills:
  - ∅3.00 & ∅3.50mm → 400 - 700 rpm.
  - ∅4.00 & ∅5.10mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.

For drilling, external cooling with saline solution is required.


frontier Ø3,30			
Descripción - Description		Ref.	
1	Fresa lanza	Lance-shaped drill	KYFoC1225
2	Fresa piloto Ø2,50 mm	Pilot drill Ø2.50 mm	KYFoC2250
3	Fresa cortical Ø3,65 mm	Cortical drill Ø3.65 mm	KYFoC4336
4	Fresa escalonada Ø2,50/Ø3,00 mm	Step drill Ø2.50/Ø3.00 mm	KYFoC5130


frontier Ø3,75			
Descripción - Description		Ref.	
1	Fresa lanza	Lance-shaped drill	KYFoC1225
2	Fresa piloto Ø2,50 mm	Pilot drill Ø2.50 mm	KYFoC2250
3	Fresa escalonada Ø2,50/Ø3,00 mm	Step drill Ø2.50/Ø3.00 mm	KYFoC5130
4	Fresa cortical Ø4,00 mm	Cortical drill Ø4.00 mm	KYFoC4340
5	Fresa escalonada Ø3,00/Ø3,50 mm	Step drill Ø3.00/Ø3.50 mm	KYFoC5135

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

gmi® frontier PEAK internal connection hexagonal implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the conical core and progressive thread profile facilitate implant insertion and act as a bone compactor allowing high primary stability to be obtained in sparse bones.

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gmi® frontier PEAK internal connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the gmi® frontier PEAK implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single hex-2.45 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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

The following is a description of the drilling sequences for the different internal connection implant models, as well as the recommended drilling conditions for their use:

- Initial drill: 1000 rpm.
- Pilot drill: 700-900 rpm.
- Step and conical drills:
  - ∅2.00 / ∅2.80 mm → 500 - 700 rpm.
  - ∅3.20 / ∅3.65 mm → 500 - 700 rpm.
  - ∅3.70 / ∅4.15 mm → 400 - 700 rpm.
  - ∅4.15 / ∅4.65 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.


For drilling, external cooling with saline solution is required.

frontier PEAK ∅3,75		
Descripción - Description		Ref.
1	Fresa inicial / Initial drill	KYFoC1401
2	Fresa piloto ∅2,50 mm / Pilot drill ∅2.50	KYFoC2250
3	Fresa escalonada ∅2,50/∅2,80 mm / Step drill ∅2.00/∅2.80 mm	KYFoC5129
4	Fresa cónica ∅3,20 mm / Conical drill ∅3.20 mm	KYFoC7032
5	Fresa cortical ∅3,65 mm / Cortical drill ∅3.65	KYFoC4337
6	Fresa escalonada ∅3,20/∅3,65mm / Step drill ∅3.20/∅3.65mm	KYFoC5137

frontier PEAK ∅4,25		
Descripción - Description		Ref.
1	Fresa inicial / Initial drill	KYFoC1401
2	Fresa piloto ∅2,50 mm / Pilot drill ∅2.50	KYFoC2250
3	Fresa escalonada ∅2,50/∅2,80 mm / Step drill ∅2.00/∅2.80 mm	KYFoC5129
4	Fresa cónica ∅3,20 mm / Conical drill ∅3.20 mm	KYFoC7032
5	Fresa cónica ∅3,70 mm / Conical drill ∅3.70 mm	KYFoC7037
6	Fresa cortical ∅4,15 mm / Cortical drill ∅4.15 mm	KYFoC4342
7	Fresa escalonada ∅3,70/∅4,15mm / Step drill ∅3.70/∅4.15mm	KYFoC5142

frontier PEAK ∅4,75		
Descripción - Description		Ref.
1	Fresa inicial / Initial drill	KYFoC1401
2	Fresa piloto ∅2,50 mm / Pilot drill ∅2.50	KYFoC2250
3	Fresa escalonada ∅2,50/∅2,80 mm / Step drill ∅2.00/∅2.80 mm	KYFoC5129
4	Fresa cónica ∅3,20 mm / Conical drill ∅3.20 mm	KYFoC7032
5	Fresa cónica ∅3,70 mm / Conical drill ∅3.70 mm	KYFoC7037
6	Fresa cónica ∅4,20 mm / Conical drill ∅4.20 mm	KYFoC7042
7	Fresa cortical ∅4,65 mm / Cortical drill ∅4.65 mm	KYFoC4347
8	Fresa escalonada ∅4,15/∅4,65mm / Step drill ∅4.15/∅4.65mm	KYFoC5147

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Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Alternative tightening with a flat tip screwdriver.
- Available in emergency heights from 2.00 to 6.00 mm.
- Apical bit to facilitate the start of threading.

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

When the use of a custom abutment is required, millable abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, millable implant abutments can be used like impression coping, threading them to the implant screw through the short impression screw an using them to make the negative of the patient's mouth for subsequent positivation on plaster.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 4.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Clinic screw ref.: KDT0C3601 included.



Plataforma Platform	Ref.
RP ●●	KDBoF3603
WP ●	KDBoF3602

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw ref.: KDT0C3601 included.
- Abutment shoulder adapted to the shape of the gum.

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Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Abutments in b/f 2.45 mm hex non-rotary version.
- Clinic screw ref: KDT0C3601 included.
- Available in three emergency heights from 1.50 to 3.50 mm.

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When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary cylinder is used, defined as technical plastic machined attachment that once attached to the crown of the implant with clinic screw, will absorb most of the masticatory stresses while avoiding overloading the implant while the osseointegration process is completed.

### Key features

- Machined in high strength technical plastic.
- Abutment shoulder adapted to the shape of the gum.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw ref.: KDT0C3601 included.

Plataforma Platform	Tipo conexión Connection type	Ref.
RP ● ●	Anti-rotatoria · Non-rotary	KDHoF3605
	Rotatoria · Rotary	KDHoF3607
WP ●	Anti-rotatoria · Non-rotary	KDHoF3606
	Rotatoria · Rotary	KDHoF3608

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
Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw ref.: KDT0C3601 included.
- Available in emergency heights of 0.60 mm, 1.50 mm and 3.00 mm.
- Non-rotary titanium base with indexation.
- Scan directly on Titanium base or alternative scan on Scan abutment.

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
\* Use the Healing abutment h=3 mm RP (KDPOC3609) and WP (KDPOC3613)

### Calcinable - Castable

Plataforma Platform	Altura total (mm) Total Height (mm)	Ref.
RP ●●	11,00	KDCoC3605
WP ●	11,00	KDCoC3606

## Base titanio gmi® frontier - gmi® frontier titanium base

### Pilar escaneado - Scanabutment

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	5,30	KDIOF3601
WP ●	5,30	KDIOF3602

### Cuerpo pilar escaneado - Scanbody

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP - WP ●●	7,50	KDIOC4003

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
Once the healing abutment has formed the mucosal tunnel and the healing abutment has been removed from this location, the corresponding ball abutment is used, which is defined as a machined attachment being fixed to the top of the implant and serving as the retention element for the removable prosthesis, with the previous addition of the retention caps that will be fixed to the top spherical part of the ball abutment.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- Pink retention cap and metal container included.
- Protection disk included.
- Tightening with b/f 2.45 mm hex wrench, part no.: KYLOC0065.
- Instruments for caps: Insertion / removal tool for caps (RH485IC).

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	1,00	KDFoF3609
	2,00	KDFoF3610
	3,00	KDFoF3611
	4,00	KDFoF3612

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.50 mm.
- One-piece straight abutments in rotary model.
- Tightened with hex keys:
  - RP: B/F 2.00 mm (KYLoC0149)
  - WP: B/F 2.50 mm (KYLoC0098)
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium as per ISO-5832 with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.


## Pilares multi-estético gmi® frontier - gmi® frontier multi-esthetic abutments


### Pilares multi-estético inclinados Multi-esthetic angled abutments

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP ● ●	17°	2,50	KDGoF3609
	17°	3,50	KDGoF3610
	30°	4,00	KDGoF3611


### Pilares multi-estético rectos Multi-esthetic straight abutments

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ● ●	1,00	KDGoC3601 ①
	1,50	KDGoC3620 ②
	2,50	KDGoC3621 ②
	3,50	KDGoC3622 ②
	4,50	KDGoC3623 ②
WP ●	1,00	KDGoC3605 ①
	1,50	KDGoC3624 ②
	2,50	KDGoC3625 ②
	3,50	KDGoC3626 ②
	4,50	KDGoC3627 ②

### Pilares cicatrización PME Healing abutments MEA


Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP ● ●	Ø 4,80	KDPoF0035
	Ø 6,00	KDPoF0037
WP ●	Ø 6,00	KDPoF0036

### Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

Plataforma Platform	Ref.
RP ● ●	KDSoF0015
WP ●	KDSoF0016

## Pilares multi-estético gmi® frontier - gmi® frontier multi-esthetic abutments

Toma impresión cubeta cerrada (CC) PME  
Closed tray (CT) impression coping MEA

Plataforma Platform	Cantidad (Ud.) Quantity (Un.)	Ref.
RP ● ●	1	KDSOF0017
	5	KDSOF5001

Cilindros provisionales PME  
Temporary cylinder MEA

Plataforma Platform	Ref.
RP ● ●	KDGoF0028
WP ●	KDGoF0029

Calcinables PME  
Castables MEA

Plataforma Platform	Ref.
RP ● ●	KDCoF0117
WP ●	KDCoF0118

Base titanio PME  
Titanium base MEA

Plataforma Platform	Ref.
RP ● ●	KDGoF0035
WP ●	KDGoF0036

Scanbody de titanio PME  
Titanium scanbody MEA

Plataforma Platform	Ref.
RP ● ●	KDIOF0002
WP ●	KDIOF0003

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
Once the insertion of the implant in the bone bed has been completed, the bioesthetic base is used, which is defined as the machined attachment that, once fixed to the coronal part of the implant, allows single restorations at the abutment level.

Placing the bioesthetic base during surgery allows the soft tissues to heal around it.

Once the tissue has healed, it serves as the basis for the rehabilitation, avoiding the connection and disconnection of the components to the implant and, therefore, keeping the biological seal intact.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- Straight two-piece bases, including the screw.
- Tightened with hex wrench: B/F 2.00 mm (KYLoC0149).
- Impression coping with screw included.
- Castable.
- Tightening torques: BE base (30 N-cm) and healing abutment (Manual: Max. 15 N-cm).

Bases bioestéticas Bioesthetic bases		
Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	1,00	KDGoF3612
	2,00	KDGoF3613
	3,00	KDGoF3614
	4,00	KDGoF3615
WP ●	1,00	KDGoF3616
	2,00	KDGoF3617
	3,00	KDGoF3618
	4,00	KDGoF3619

Pilar cicatrización BE BE healing abutment		
Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP - WP ●●	Ø 4,80	KDPoC0038

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


Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

### Key features


- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.


Réplicas de implante Implant replicas	
Plataforma Platform	Ref.
RP ●●	KDRoC3601
WP ●	KDRoC3602

Réplicas digitales Digital replicas	
Plataforma Platform	Ref.
RP ●●	KDRoF3603
WP ●	KDRoF3604

Réplicas para pilares multi-estético MEA replicas	
Plataforma Platform	Ref.
RP ●●	KDRoC0012
WP ●	KDRoC0013

Réplicas digitales para pilares multi-estético MEA digital replicas	
Plataforma Platform	Ref.
RP ●●	KDRoF0014
WP ●	KDRoF0015

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Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, are used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

### Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- CT impression coping whole formed by CT impression coping + hexagonal abutment + screw.

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

To replace a dental crown with a restoration dental prosthesis, a castable is used, which is defined as a meltable plastic machined attachment that, once positioned over the coronal part of the internal connection implant, will be the initial basis to model the form of the final structure, to be cast in the corresponding metal.

### Key features

- Machined castables in white meltable plastic material.
- Physiologically inert and excellent dimensional stability.
- Available in rotary and non-rotary version.
- Clinic screw ref.: KDTC3601 included.

Plataforma Platform	Tipo de conexión Connection type	Ref.
RP ● ●	Anti-rotatoria • Non-rotary	KDC0F3611
	Rotatoria • Rotary	KDC0F3612
WP ●	Anti-rotatoria • Non-rotary	KDC0F3613
	Rotatoria • Rotary	KDC0F3614

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
The CrCo bases consist of a machined CrCo interface and a castable chimney. This configuration allows for greater dimensional stability at the contact surfaces of the implant-abutment connection.

They are indicated for single or multiple screw-retained restorations either by conventional cast-on technique or by laser welding.

Plataforma Platform	Tipo de conexión Connection type	Ref.
RP ●●	Rotatoria • Rotary	TS30017
	Anti-rotatoria • Non-rotary	TS31017
WP ●	Rotatoria • Rotary	TS30018
	Anti-rotatoria • Non-rotary	TS31018

## Tornillos clínica y laboratorio gmi® frontier - gmi® frontier clinic and laboratory screws

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
Clinic screws are defined as machined devices used for the final anchoring of the abutments, the millable implant abutment or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of castables to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo clínica Clinic screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ●●	○ 1,20 mm	KDT0C3601
	⊗ TX	KDT0C3609
Tornillo laboratorio Laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ●●	○ 1,20 mm	KDT0C3602
	⊗ TX	KDT0C3610

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MEA and BE clinical screws are defined as machined devices used for the final anchoring of the ME and BE abutments, while laboratory screws are the machined devices used for the temporary anchoring of ME and BE abutments to the replica during the fabrication of the prostheses in the laboratory.

### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm and 2.00 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo PME MEA screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP ●●	⊕ 1,20 mm	KDToC0056
	⊕ TX	KDToC0070
WP ●	⊕ 1,20 mm	KDToC0060
	⊕ TX	KDToC0071

Tornillo laboratorio PME MEA laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP ●●	⊕ 1,20 mm	KDToC0057
	⊕ TX	KDToC0072
WP ●	⊕ 1,20 mm	KDToC0058
	⊕ TX	KDToC0073

Tornillo BE BE screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ●●	⊕ 1,20 mm	KDToC4016
	⊕ TX	KDToC0074

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gmi® frontier/avantgard surgical kit (Ref. KYCoF3053) consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of gmi frontier and avantgard dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi® surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134°C.

To complement this kit, gmi® has the PEAK drill Kit (Ref. KYCoF3056 - see instrumental catalogue), allowing it to be easily expanded for the placement of 2 additional implant systems (frontier PEAK and avantgard PEAK).

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Ref. KYCoF3058 - Kit for gmi® frontier system only, without the drill and wrenches exclusive to the avantgard system (see table, highlighted in blue).

Descripción - Description		Un.	Ref.
Fresa lanza	Lance-shaped drill	1	KYFoC1225
Fresa piloto ø2,50 mm	Pilot drill ø2.50 mm	1	KYFoC2250
Fresa escalonada ø2,50 / ø3,00 mm	Step drill ø2.50 / ø3.00 mm	1	KYFoC5130
Fresa escalonada ø3,00 / ø3,50 mm	Step drill ø3.00 / ø3.50 mm	1	KYFoC5135
Fresa escalonada ø3,50 / ø4,00 mm	Step drill ø3.50 / ø4.00 mm	1	KYFoC5140
Fresa escalonada ø4,00 / ø4,50 mm	Step drill ø4.00 / ø4.50 mm	1	KYFoC5145
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4336
Fresa cortical ø3,20 mm	Cortical drill ø3.20 mm	1	KYFoC4332
Fresa cortical ø4,00 mm	Cortical drill ø4.00 mm	1	KYFoC4340
Fresa cortical ø4,50 mm	Cortical drill ø4.50 mm	1	KYFoC4345
Fresa cortical ø4,90 mm	Cortical drill ø4.90 mm	1	KYFoC4349
Fresa cortical ø5,90 mm	Cortical drill ø5.90 mm	1	KYFoC4359
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave-transportador HEX-2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYLoF0096
Llave-transportador HEX-2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYLoF0111
Llave-transportador HEX-2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYLoF0091
Llave Avantgard NP/RP IP corta	Short Avantgard NP/RP TI wrench	1	KYLoF0157
Llave Avantgard NP/RP IP larga	Long Avantgard NP/RP TI wrench	1	KYLoF0162
Llave Avantgard WP IP corta	Short Avantgard WP TI wrench	1	KYLoF0168
Llave Avantgard WP IP larga	Long Avantgard WP TI wrench	1	KYLoF0169
Llave Avantgard NP/RP CA corta	Short Avantgard NP/RP HP wrench	1	KYLoF0151
Llave Avantgard WP CA corta	Short Avantgard WP HP wrench	1	KYLoF0166
Llave carraca corta HEX-1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Llave carraca larga HEX-1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLoF0129
Alargador fresas	Drill extender	1	KYLoC0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLoC0078
Indicador profundidad ø2,50 mm	Depth gauge ø2.50 mm	1	KYLoC0165

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gmi® frontier guided surgery kit (Ref. KYCoF3047) consists of an autoclavable case made of high strength resin, containing all the necessary material for performing totally guided surgeries with gmi® frontier implants of  $\varnothing 3.30$  /  $\varnothing 3.75$   $\varnothing 4.25$  and  $\varnothing 4.75$  mm and lengths between 8 and 15 mm.

To complement this kit, gmi® has the PEAK guided surgery drill Kit (Ref. KYCoF3052 - see instrumental catalogue), allowing it to be easily expanded for the placement of 1 additional implant system (frontier PEAK).

Descripción - Description		Un.	Ref.
Fresa Pin $\varnothing 1,50$ mm	Pin drill $\varnothing 1.50$ mm	1	KYFoC6060
Trepante gingival CA $\varnothing 5,10$ mm	HP gingival punch $\varnothing 5.10$ mm	1	KYFoC6066
Fresa inicial $\varnothing 2,5$ mm	Initial drill $\varnothing 2.5$ mm	1	KYFoC6101
Fresa piloto DLC $\varnothing 2,5$ L8 mm	DLC pilot drill $\varnothing 2.5$ L8 mm	1	KYFoC6201
Fresa piloto DLC $\varnothing 2,5$ L10 mm	DLC pilot drill $\varnothing 2.5$ L10 mm	1	KYFoC6202
Fresa piloto DLC $\varnothing 2,5$ L11,5 mm	DLC pilot drill $\varnothing 2.5$ L11.5 mm	1	KYFoC6203
Fresa piloto DLC $\varnothing 2,5$ L13 mm	DLC pilot drill $\varnothing 2.5$ L13 mm	1	KYFoC6204
Fresa piloto DLC $\varnothing 2,5$ L15 mm	DLC pilot drill $\varnothing 2.5$ L15 mm	1	KYFoC6205
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L8 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L8 mm	1	KYFoC6506
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L10 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L10 mm	1	KYFoC6507
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L11,5 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L11.5 mm	1	KYFoC6508
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L13 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L13 mm	1	KYFoC6509
Fresa escalonada DLC $\varnothing 2,5/\varnothing 3$ L15 mm	DLC step drill $\varnothing 2.5/\varnothing 3$ L15 mm	1	KYFoC6510
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L8 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L8 mm	1	KYFoC6511
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L10 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L10 mm	1	KYFoC6512
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L11,5 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L11.5 mm	1	KYFoC6513
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L13 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L13 mm	1	KYFoC6514
Fresa escalonada DLC $\varnothing 3/\varnothing 3,5$ L15 mm	DLC step drill $\varnothing 3/\varnothing 3.5$ L15 mm	1	KYFoC6515
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L8 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L8 mm	1	KYFoC6516
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L10 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L10 mm	1	KYFoC6517
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L11,5 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L11.5 mm	1	KYFoC6518
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L13 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L13 mm	1	KYFoC6519
Fresa escalonada DLC $\varnothing 3,5/\varnothing 4$ L15 mm	DLC step drill $\varnothing 3.5/\varnothing 4$ L15 mm	1	KYFoC6520
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L8 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L8 mm	1	KYFoC6521

Descripción - Description		Un.	Ref.
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L10 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L10 mm	1	KYFoC6522
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L11,5 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L11.5 mm	1	KYFoC6523
Fresa escalonada DLC $\varnothing 4/\varnothing 4,5$ L13 mm	DLC step drill $\varnothing 4/\varnothing 4.5$ L13 mm	1	KYFoC6524
Fresa cortical DLC $\varnothing 3,2$ mm	DLC cortical drill $\varnothing 3.2$ mm	1	KYFoC6401
Fresa cortical DLC $\varnothing 3,55$ mm	DLC cortical drill $\varnothing 3.55$ mm	1	KYFoC6402
Fresa cortical DLC $\varnothing 3,9$ mm	DLC cortical drill $\varnothing 3.9$ mm	1	KYFoC6403
Fresa cortical DLC $\varnothing 4,4$ mm	DLC cortical drill $\varnothing 4.4$ mm	1	KYFoC6404
Fresa cortical DLC $\varnothing 4,9$ mm	DLC cortical drill $\varnothing 4.9$ mm	1	KYFoC6405
Llave carraca corta HEX-1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Adaptador HEX-4,00 mm LLCIP	HEX-4.00 mm adapter TIRW	1	KYLoC6054
Adaptador HEX-4,00 mm CA	HEX-4.00 mm adapter HP	1	KYLoC6055
Extractor transportadores	Mounters extraction wrench	1	KYLoC6115
Transportadores RP/WP + tornillo	Implant carriers RP/WP + screw	4	KYLoF6103
Llave directa RP/WP CA	Direct HP wrench RP/WP	1	KYLoF6105
Llave directa RP/WP LLCIP	Direct TIRW wrench RP/WP	1	KYLoF6106
Pin lateral férula	Lateral pin	3	KYLoC6059
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave transportador plana 45°	45° implant carrier manual wrench	1	KYLoC0093

\*Todas las referencias se venden por separado como recambios  
\*All references are sold separately as spare parts

Accesorios - Accessories		Un.	Ref.
Casquillo pin lateral férula	Sleeve for lateral pin	5	KYLoF6061
Casquillo férula $\varnothing 5,10$ mm	Sleeve $\varnothing 5.10$ mm	5	KYLoF6102

## Kit quirúrgico gmi® PEAK - gmi® PEAK surgical kit

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gmi® PEAK surgical kit (Ref. KYCoF3054) consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of frontier PEAK and avantgard PEAK dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi® surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134°C.

To complement this kit, gmi® has the **frontier/avantgard drill Kit** (Ref. KYCoF3055 - *see instrumental catalogue*), allowing it to be easily expanded for the placement of 2 additional implant systems (frontier and avantgard).

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Ref. KYCoF3060 - Kit for gmi® frontier PEAK system only, without the drill and wrenches exclusive to the avantgard PEAK system (*see table, highlighted in blue*).

Descripción - Description		Un.	Ref.
Fresa inicial ø2,0 mm	Initial drill ø2.0 mm	1	KYFoC1401
Fresa piloto ø2,50 mm	Pilot drill ø2.50 mm	1	KYFoC2250
Fresa cónica ø3,20 mm	Conical drill ø3.20 mm	1	KYFoC7032
Fresa cónica ø3,70 mm	Conical drill ø3.70 mm	1	KYFoC7037
Fresa cónica ø4,20 mm	Conical drill ø4.20 mm	1	KYFoC7042
Fresa cortical ø3,20 mm	Cortical drill ø3.20 mm	1	KYFoC4333
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4337
Fresa cortical ø4,15 mm	Cortical drill ø4.15 mm	1	KYFoC4342
Fresa cortical ø4,65 mm	Cortical drill ø4.65 mm	1	KYFoC4347
Fresa escalonada ø2,50 / ø2,80 mm	Step drill ø2.50 / ø2.80 mm	1	KYFoC5129
Fresa escalonada ø3,20 / ø3,65 mm	Step drill ø3.20 / ø3.65 mm	1	KYFoC5137
Fresa escalonada ø3,70 / ø4,15 mm	Step drill ø3.70 / ø4.15 mm	1	KYFoC5142
Fresa escalonada ø4,15 / ø4,65 mm	Step drill ø4.15 / ø4.65 mm	1	KYFoC5147
Indicador profundidad ø2,50 mm	Depth gauge ø2.50 mm	1	KYLoC0165
Indicador profundidad cónico ø3,20 mm	Conical depth gauge ø3.20 mm	1	KYLoC0159
Indicador profundidad cónico ø3,70 mm	Conical depth gauge ø3.70 mm	1	KYLoC0160
Indicador profundidad cónico ø4,20 mm	Conical depth gauge ø4.20 mm	1	KYLoC0161
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave-transportador HEX - 2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYLoF0091
Llave-transportador HEX - 2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYLoF0096
Llave-transportador HEX - 2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYLoF0111
Llave Avantgard NP/RP CA corta	Short Avantgard NP/RP HP wrench	1	KYLoF0151
Llave Avantgard NP/RP IP corta	Short Avantgard NP/RP TI wrench	1	KYLoF0157
Llave Avantgard NP/RP IP larga	Long Avantgard NP/RP TI wrench	1	KYLoF0162
Alargador fresas	Drill extender	1	KYLoC0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLoC0078
Llave carraca corta HEX-1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Llave carraca larga HEX-1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLoF0129

## Kit cirugía guiada gmi® frontier PEAK - gmi® frontier PEAK guided surgery kit

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gmi® frontier PEAK guided surgery kit (Ref. KYCoF3048) consists of an autoclavable case made of high strength resin, containing all the necessary material for performing totally guided surgeries with gmi® frontier PEAK implants of  $\varnothing 3.75$   $\varnothing 4.25$  and  $\varnothing 4.75$  mm and lengths between 8 and 15 mm.

To complement this kit, gmi® has the frontier/avantgard guided surgery drill Kit (Ref. KYCoF3051 - *see instrumental catalogue*), allowing it to be easily expanded for the placement of 1 additional implant system (frontier).

Descripción - Description		Un.	Ref.
Fresa Pin $\varnothing 1,50$ mm	Pin drill $\varnothing 1.50$ mm	1	KYFoC6060
Trepante gingival CA $\varnothing 5,10$ mm	HP gingival punch $\varnothing 5.10$ mm	1	KYFoC6066
Fresa inicial $\varnothing 2,5$ mm	Initial drill $\varnothing 2.5$ mm	1	KYFoC6101
Fresa piloto DLC $\varnothing 2,5$ L8 mm	DLC pilot drill $\varnothing 2.5$ L8 mm	1	KYFoC6206
Fresa piloto DLC $\varnothing 2,5$ L10 mm	DLC pilot drill $\varnothing 2.5$ L10 mm	1	KYFoC6207
Fresa piloto DLC $\varnothing 2,5$ L11,5 mm	DLC pilot drill $\varnothing 2.5$ L11.5 mm	1	KYFoC6208
Fresa piloto DLC $\varnothing 2,5$ L13 mm	DLC pilot drill $\varnothing 2.5$ L13 mm	1	KYFoC6209
Fresa piloto DLC $\varnothing 2,5$ L15 mm	DLC pilot drill $\varnothing 2.5$ L15 mm	1	KYFoC6210
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L8 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L8 mm	1	KYFoC6525
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L10 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L10 mm	1	KYFoC6526
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L11,5 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L11.5 mm	1	KYFoC6527
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L13 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L13 mm	1	KYFoC6528
Fresa escalonada DLC $\varnothing 2,5/\varnothing 2,8$ L15 mm	DLC step drill $\varnothing 2.5/\varnothing 2.8$ L15 mm	1	KYFoC6529
Fresa cónica DLC $\varnothing 3,2$ L8 mm	DLC conical drill $\varnothing 3.2$ L8 mm	1	KYFoC6706
Fresa cónica DLC $\varnothing 3,2$ L10 mm	DLC conical drill $\varnothing 3.2$ L10 mm	1	KYFoC6707
Fresa cónica DLC $\varnothing 3,2$ L11,5 mm	DLC conical drill $\varnothing 3.2$ L11.5 mm	1	KYFoC6708
Fresa cónica DLC $\varnothing 3,2$ L13 mm	DLC conical drill $\varnothing 3.2$ L13 mm	1	KYFoC6709
Fresa cónica DLC $\varnothing 3,2$ L15 mm	DLC conical drill $\varnothing 3.2$ L15 mm	1	KYFoC6710
Fresa cónica DLC $\varnothing 3,7$ L8 mm	DLC conical drill $\varnothing 3.7$ L8 mm	1	KYFoC6711
Fresa cónica DLC $\varnothing 3,7$ L10 mm	DLC conical drill $\varnothing 3.7$ L10 mm	1	KYFoC6712
Fresa cónica DLC $\varnothing 3,7$ L11,5 mm	DLC conical drill $\varnothing 3.7$ L11.5 mm	1	KYFoC6713
Fresa cónica DLC $\varnothing 3,7$ L13 mm	DLC conical drill $\varnothing 3.7$ L13 mm	1	KYFoC6714
Fresa cónica DLC $\varnothing 3,7$ L15 mm	DLC conical drill $\varnothing 3.7$ L15 mm	1	KYFoC6715
Fresa cónica DLC $\varnothing 4,2$ L8 mm	DLC conical drill $\varnothing 4.2$ L8 mm	1	KYFoC6716

Descripción - Description		Un.	Ref.
Fresa cónica DLC $\varnothing 4,2$ L10 mm	DLC conical drill $\varnothing 4.2$ L10 mm	1	KYFoC6717
Fresa cónica DLC $\varnothing 4,2$ L11,5 mm	DLC conical drill $\varnothing 4.2$ L11.5 mm	1	KYFoC6718
Fresa cónica DLC $\varnothing 4,2$ L13 mm	DLC conical drill $\varnothing 4.2$ L13 mm	1	KYFoC6719
Fresa cortical DLC $\varnothing 3,2$ mm	DLC cortical drill $\varnothing 3.2$ mm	1	KYFoC6406
Fresa cortical DLC $\varnothing 3,65$ mm	DLC cortical drill $\varnothing 3.65$ mm	1	KYFoC6407
Fresa cortical DLC $\varnothing 4,15$ mm	DLC cortical drill $\varnothing 4.15$ mm	1	KYFoC6408
Fresa cortical DLC $\varnothing 4,65$ mm	DLC cortical drill $\varnothing 4.65$ mm	1	KYFoC6409
Llave carraca corta Hex-1,20 mm	Short Hex-1.20 mm ratchet wrench	1	KYLoF0128
Adaptador Hex-4,00 mm LLCIP	Hex-4.00 mm adapter TIRW	1	KYLoC6054
Adaptador Hex-4,00 mm CA	Hex-4.00 mm adapter HP	1	KYLoC6055
Extractor transportadores	Mounters extraction wrench	1	KYLoC6115
Transportadores RP/WP + tornillo	Implant carriers RP/WP + screw	4	KYLoF6103
Llave directa RP/WP CA	Direct HP wrench RP/WP	1	KYLoF6105
Llave directa RP/WP LLCIP	Direct TIRW wrench RP/WP	1	KYLoF6106
Pin lateral férula	Lateral pin	3	KYLoC6059
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave transportador plana 45°	45° implant carrier manual wrench	1	KYLoC0093

\*Todas las referencias se venden por separado como recambios  
 \*All references are sold separately as spare parts

Accesorios - Accessories		Un.	Ref.
Casquillo pin lateral férula	Sleeve for lateral pin	5	KYLoF6061
Casquillo férula $\varnothing 5,10$ mm	Sleeve $\varnothing 5.10$ mm	5	KYLoF6102

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