SISTEMA DE IMPLANTES  
IMPLANT SYSTEM

[www.gmidental.com](http://www.gmidental.com)

Implantes dentales de conexión externa  
External connection dental implants

Made in Spain

## Índice - Summary

### CARACTERÍSTICAS GENERALES

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### KIT QUIRÚRGICO GMI PHOENIX

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### KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

### IMPLANT SYSTEM

- **gmi phoenix** dental implants 8
- Insertion sequence for **gmi phoenix** implants 9
- **gmi phoenix** implants range 10
- Drilling sequences for **gmi phoenix** implants 11

### ATTACHMENTS

- **gmi phoenix** healing abutments 13
- **gmi phoenix** esthetic angled abutments 15
- **gmi phoenix** straight abutments 16
- **gmi phoenix** implant carrier abutments 17
- **gmi phoenix** ball abutments 18
- **gmi phoenix** multi-esthetic abutments 20
- **gmi phoenix** replicas 23
- **gmi phoenix** impression coping 24
- **gmi phoenix** castables 26
- **gmi phoenix** clinic and laboratory screws 28
- **gmi phoenix** MEA screws 29
- Recommended tightening torques 30

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


To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.





SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: ultrastructural details of cells in biofilm on the implant.

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The **gmi phoenix** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

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
The **gmi phoenix** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi phoenix** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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

To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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
In the chart are detailed all the symbols appearing on the **gmi phoenix** implant system's labelling and packaging and their corresponding description.

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\*Minimum requirements for downloading and opening the electronic instructions for use:  
- Stable internet connectivity.  
- A device with a modern operating system that supports PDF viewing.  
- Adequate storage space on your device.  
- A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.  
- A PDF reader application installed on your device.

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
gmi phoenix external connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

### Características generales

#### Key features

#### ① Rosca interior

Estrecha: M1,8 x 0,35 mm  
 Estándar: M2,0 x 0,40 mm  
 Ancha: M2,0 x 0,40 y M2,5 x 0,45 mm  
 Extra-Ancha: M2,5 x 0,45 mm

#### ② Hexágono externo

Estrecha: E/C 2,40 mm x 1 mm  
 Estándar: E/C 2,70 mm x 0,7 mm  
 Ancha: E/C 2,70 mm x 0,7 mm  
           E/C 3,40 mm x 0,7 mm  
 Extra-Ancha: E/C 3,40 x 0,9 mm

#### ③ Plataforma

#### ④ Rosca Exterior

Paso 1,60 mm  
 Tratamiento superficial ADS

#### ⑤ Fresados auto-roscantes

Longitud Variable  
 3 fresados a 120°

#### ⑥ Longitud

#### ① Internal thread

Narrow: M1.8 x 0.35 mm  
 Standard: M2.0 x 0.40 mm  
 Wide: M2.0 x 0.40 & M2.5 x 0.45 mm  
 Extra wide: M2.5 x 0.45 mm

#### ② External hexagon

Narrow: B/F 2.40 mm x 1 mm  
 Standard: B/F 2.70 mm x 0.7 mm  
 Wide: B/F 2.70 mm x 0.7 mm  
           B/F 3.40 mm x 0.7 mm  
 Extra wide: B/F 3.40 x 0.9 mm

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**gmi phoenix** external connection implants are packed in a screen-printed cardboard box, with 3 traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi phoenix** implants consists of a titanium support which only houses the implant-implant carrier set, prepared to be inserted directly into the mouth, using a single hex-3.00 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.\*

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
\* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the tab on the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the HEX-3.00 mm implant carrier wrench until a slight retention is noticed, addressing the notches on the key with the faces of the implant hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the crestal or subcrestal level. Remove the key.
9. Unscrew the clinic screw and dismount the implant carrier by turning slightly if necessary.
10. Remove the cover screw from the primary packaging plug, using the HEX-0.80 mm ratchet wrench for Ø3.30 platform and HEX-1.20 mm for the other platforms.
11. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove key and suture incision.

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


Following is a description of the drilling sequences for the different external connection implant models, as well as the recommended drilling conditions for their use.

- Lance-shaped: 1200 - 1500 rpm
- Pilot drill: 900 - 1200 rpm
- Guide drill: 800 rpm
- Twist drills:
  - $\varnothing 2.80$  mm  $\longrightarrow$  500 - 700 rpm.
  - $\varnothing 3.00$  and  $\varnothing 3.50$  mm  $\longrightarrow$  400 - 700 rpm.
  - $\varnothing 4.25$  mm  $\longrightarrow$  400 - 600 rpm.
  - $\varnothing 5.10$  and  $\varnothing 5.40$  mm  $\longrightarrow$  300 - 500 rpm.
- For drilling, external cooling with saline solution is required.


phoenix ø3,30		
	Descripción - Description	Ref.
1	Fresa lanceolada Lance-shaped drill	KYFoC1225
2	Fresa piloto Pilot drill	KYFoC2221
3	Fresa guía Guide drill	KYFoC3227
4	Fresa $\varnothing 2,80$ mm Twist drill $\varnothing 2.80$ mm	KYFoC0128


phoenix ø3,75		
	Descripción - Description	Ref.
1	Fresa lanceolada Lance-shaped drill	KYFoC1225
2	Fresa piloto Pilot drill	KYFoC2221
3	Fresa guía Guide drill	KYFoC3227
4	Fresa $\varnothing 2,80$ mm Twist drill $\varnothing 2.80$ mm	KYFoC0128
5	Fresa $\varnothing 3,25$ mm Twist drill $\varnothing 3.25$ mm	KYFoC0132

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
Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Alternative tightening with a flat tip screwdriver.
- Available in emergency heights from 2.00 to 5.00 mm.
- Slightly tapered in the head to facilitate extraction.
- Apical bit to facilitate the start of threading.

### Pilares de cicatrización para plataforma $\varnothing 3,30$ mm Healing abutments for $\varnothing 3.30$ mm platform

Hex. impl.	Altura (mm) Height (mm)	Ref.
$\varnothing 2,40$ mm	2,00	KDPOC0023
	3,00	KDPOC0024
	4,00	KDPOC0025
	5,00	KDPOC0026

### Pilares de cicatrización para plataforma $\varnothing 4,10$ mm Healing abutments for $\varnothing 4.10$ mm platform

Hex. impl.	Altura (mm) Height (mm)	Ref.
$\varnothing 2,70$ mm	2,00	KDPOC0005
	3,00	KDPOC0006
	4,00	KDPOC0007
	5,00	KDPOC0008

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

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
### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- 12-position interior for increased flexibility.
- Clinic screw included.
- Abutment shoulder adapted to the shape of the gum.

### Pilares inclinados estéticos para plataforma ø3,30 mm Esthetic angled abutments for ø3.30 mm platform

Hex. impl.	Tipo conexión Connection type	Ang. (°)	Ref.
○ 2,40 mm	Antir. • Non-rot.	15°	KDMoF0060
	Antir. • Non-rot.	20°	KDMoF0069

### Pilares inclinados estéticos para plataforma ø4,10 mm Esthetic angled abutments for ø4.10 mm platform

Hex. impl.	Tipo conexión Connection type	Ang. (°)	Ref.
○ 2,70 mm	Antir. • Non-rot.	15°	KDMoF0002
	Antir. • Non-rot.	20°	KDMoF0070

### Pilares inclinados estéticos para plataforma ø5,10 mm Esthetic angled abutments for ø5.10 mm platform

Hex. impl.	Tipo conexión Connection type	Ang. (°)	Ref.
○ 2,70 mm	Antir. • Non-rot.	15°	KDMoF0008
	Antir. • Non-rot.	20°	KDMoF0071
○ 3,40 mm	Antir. • Non-rot.	15°	KDMoF0014
	Antir. • Non-rot.	20°	KDMoF0072

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Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Straight abutments in hex non-rotary version.
- Clinic screw included.

### Pilares rectos para plataforma ø3,30 mm Straight abutments for ø3.30 mm platform

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When the use of a custom abutment is required, carrier abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, millable transfer abutments can be used like impression coping, threading them to the implant screw through the short impression screw and using them to make the negative of the patient's mouth for subsequent positivation on plaster.

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### Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 4.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Clinic screw included.

### Pilares transportador para plataformas $\varnothing 3,30 / 4,10 / 5,10 / 6,00$ mm Implant carrier abutments for $\varnothing 3.30 / 4.10 / 5.10 / 6.00$ mm platforms

Hex. impl.	Plataforma Platform	Ref.
$\varnothing 2,40$ mm	$\varnothing 3.30$	KDBoF0004
$\varnothing 2,70$ mm	$\varnothing 4.10$	KDBoF0001
	$\varnothing 5.10$	KDBoF0002
$\varnothing 3,40$ mm	$\varnothing 5.10$	KDBoF0003
	$\varnothing 6.00$	KDBoF0005

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Once the healing abutment has formed the mucosal tunnel and the healing abutment has been removed from this location, the corresponding ball abutment is used, which is defined as a machined attachment being fixed to the top of the implant and serving as the retention element for the removable prosthesis, with the previous addition of the retention caps that will be fixed to the top spherical part of the ball abutment.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.50 to 3.00 mm.
- Pink retention cap and metal container included.
- Protection disk included.
- Tightening with b/f 2.45 mm hex wrench, part no.: KYLOCO065.
- Instruments for caps: Insertion / removal tool for caps (RH485IC).

### Pilares bola para plataforma $\varnothing 4,10$ mm Ball abutment for $\varnothing 4.10$ mm platform

Hex. impl.	Altura (mm) Height (mm)	Ref.
$\varnothing 2,70$ mm	3,00	KDFoF0045
	2,50	KDFoF0046
	2,00	KDFoF0047
	1,50	KDFoF0048

### Pilares bola para plataforma $\varnothing 5,10$ mm Ball abutment for $\varnothing 5.10$ mm platform

Hex. impl.	Altura (mm) Height (mm)	Ref.
$\varnothing 3,40$ mm	3,00	KDFoF0049
	2,50	KDFoF0050
	2,00	KDFoF0051
	1,50	KDFoF0052
$\varnothing 2,70$ mm	3,00	KDFoF0053
	2,50	KDFoF0054
	2,00	KDFoF0055
	1,50	KDFoF0056

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.50 to 4.50 mm.
- One-piece straight abutments in rotary model.
- Tightened with hex keys:  
 PL-Ø3.30 and Ø4.10: B/F 2.00 mm (KYLOCo149)  
 PL-Ø5.10: B/F 2.50 mm (KYLOCo098)
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.

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


Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

### Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.

### Réplicas de implante para plataformas $\varnothing 3,30$ / $4,10$ / $5,10$ / $6,00$ mm Implant replicas for $\varnothing 3.30$ / $4.10$ / $5.10$ / $6.00$ mm platforms

Hex. impl.	Plataforma (mm) Platform (mm)	Ref.
$\varnothing 2,40$ mm	$\varnothing 3,30$	KDRoC0010
$\varnothing 2,70$ mm	$\varnothing 4,10$	KDRoC0002
	$\varnothing 5,10$	KDRoC0006
$\varnothing 3,40$ mm	$\varnothing 5,10$	KDRoC0008
	$\varnothing 6,00$	KDRoC0011

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Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, are used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

### Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.

### Tomas impresión cubeta abierta (CA) plataforma Ø3,30 mm Open tray (OT) impression coping for Ø3.30 mm platform

Hex. impl.	Tornillo Screw	Ref.
○ 2,40 mm	Largo • Long	KDS0F0007
	Corto • Short	KDS0F0012

### Tomas impresión cubeta abierta (CA) plataforma Ø4,10 mm Open tray (OT) impression coping for Ø4.10 mm platform

Hex. impl.	Tornillo Screw	Ref.
○ 2,70 mm	Largo • Long	KDS0F0002
	Corto • Short	KDS0F0009

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**Tomas impresión cubeta cerrada (CC)**  
 Closed tray (CT) impression coping

Cantidad (Ud.) Quantity (Un.)	Ref.
5	KDS0F0014

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To replace a dental crown with a restoration dental prosthesis, a castable is used, which is defined as a meltable plastic machined attachment that, once positioned over the coronal part of the external connection implant, will be the initial basis to model the form of the final structure, to be cast in the corresponding metal.

### Key features

- Machined castables in white meltable plastic material.
- Physiologically inert and excellent dimensional stability.
- Available in rotary and non-rotary version.

### Calcinables para plataforma $\varnothing 3,30$ mm Castables for $\varnothing 3.30$ mm platform

Hex. Impl.	Tipo conexión Connection type	Ref.
$\varnothing 2,40$ mm	Anti-rotatoria • Non-rotary	KDCoF0115
	Rotatoria • Rotary	KDCoF0116

### Calcinables para plataforma $\varnothing 4,10$ mm Castables for $\varnothing 4.10$ mm platform

Hex. Impl.	Tipo conexión Connection type	Ref.
$\varnothing 2,70$ mm	Anti-rotatoria • Non-rotary	KDCoF0101
	Rotatoria • Rotary	KDCoF0105

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Clinic screws are defined as machined devices used for the final anchoring of the abutments, the transfer or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of castables to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

### Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading.

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MEA screws are defined as machined devices used for the final anchoring of the multi-esthetic abutments.

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- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.

Tornillos PME MEA screws	
Plataforma Platform	Ref.
ø4,80	KDToC0056
ø6,00	KDToC0060

Tornillos laboratorio PME MEA laboratory screws	
Plataforma Platform	Ref.
ø4,80	KDToC0057
ø6,00	KDToC0058

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**gmi phoenix surgical kit** consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi phoenix** external connection dental implant system. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

**gmi surgical kits** are distributed in cartons, made of high strength technical polymer, and should be steam autoclaved before use at a maximum temperature of 134°C.

	Descr.	Un.	Ref.
Fresa lanceolada	Lance-shaped drill	1	KYFoC1225
Fresa piloto ø2,00 mm	Pilot drill ø2.00 mm	1	KYFoC2221
Fresa guía ø2,00 / ø2,70 mm	Guide drill ø2.00 / ø2.70 mm	1	KYFoC3227
Fresa final ø2,80 mm	Twist drill ø2.80 mm	1	KYFoC0128
Fresa final ø3,00 mm	Twist drill ø3.00 mm	1	KYFoC0130
Fresa final ø3,25 mm	Twist drill ø3.25 mm	1	KYFoC0132
Fresa final ø3,50 mm	Twist drill ø3.50 mm	1	KYFoC0135
Fresa final ø4,25 mm	Twist drill ø4.25 mm	1	KYFoC0143
Fresa final ø5,10 mm	Twist drill ø5.10 mm	1	KYFoC0151
Fresa final ø5,40 mm	Twist drill ø5.40 mm	1	KYFoC0154
Llave transportador plana 45°	Manual wrench 45°	1	KYLoC0093
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave-transportador HEX - 3,00 mm IP corta	Short HEX-3.00 mm implant carrier TI wrench	1	KYLoF0123
Llave-transportador HEX - 3,00 mm IP larga	Long HEX-3.00 mm implant carrier TI wrench	1	KYLoF0124
Llave-transportador HEX - 3,00 mm CA	HEX-3.00 mm implant carrier HP wrench	1	KYLoF0122
Llave fija HEX - 1,20 mm manual	HEX-1.20 mm fixed manual wrench	1	KYLoC0011
Llave carraca corta HEX - 1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Llave carraca larga HEX - 1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLoF0129
Llave fija HEX - 0,80 mm manual	HEX-0.80 mm fixed manual wrench	1	KYLoC0092
Alargador fresas	Drill extender	1	KYLoC0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLoC0078
Indicador profundidad ø2,80 mm	Depth gauge ø2.80 mm	1	KYLoC0094
Indicador profundidad ø3,50 mm	Depth gauge ø3.50 mm	1	KYLoC0095

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