

GMI STUDY. COMPARATIVE PROTEOMIC STUDY BETWEEN **GMI** and **STRAUMANN** DENTAL IMPLANTS




Hypothesis. Proteomic study of first protein layer

Post-implantation, the biomaterial becomes in contact with the blood,




Protein adsorption on its surface

(Andersson, 2005).



Protein importance on biomaterial





The first layer of adsorbed proteins will define



The biological response to the material


(Engberg, 2015)



Cell-Protein-Material interaction

Regeneration

GMI proteomic comparative study



Proteomics: methodology

1



Incubation on human serum (3h 37°C)

Consecutive washes to eliminate the non-adsorbed proteins

H₂O Mili-Q
100mM NaCl Tris-HCl pH 7.1
N=4 (4 implants)

Biomaterial surface-attached protein elution

4% SDS, 100 mM DTT in 0.5M TEAB

2



Final elution

3

Quantification and identification:
Mass Spectrometry (MS/MS)

Data analysis: PROGENESIS

DAVID (Database for Annotation, Visualization and Integrated Discovery)

Biological data base of protein families, useful to protein function identification and classification

4



DAVID Bioinformatics Resources 6.7
National Institute of Allergy and Infectious Diseases (NIAID), NIH

COMPARATIVE PROTEOMIC STUDY:

ST



GMI



Comparison GMI vs ST

Proteins with higher affinity with GMI implants



			GMI/ST	
Accession	Peptide count	Description	p value	Ratio
PROC_HUMAN	8	Vitamin K-dependent protein C	1,17E-06	27,21
FA10_HUMAN	8	Coagulation factor X	1,86E-04	17,21
FA7_HUMAN	3	Coagulation factor VII	1,04E-04	16,27
TAGL2_HUMAN	2	Transgelin-2	2,40E-05	11,54
SAA1_HUMAN	2	Serum amyloid A-1 protein	1,46E-02	5,48
FA9_HUMAN	7	Coagulation factor IX	2,44E-04	5,47
APOC1_HUMAN	3	Apolipoprotein C-I	3,58E-04	5,05
FCN3_HUMAN	2	Ficolin-3	5,42E-03	4,65
CAMP_HUMAN	2	Cathelicidin antimicrobial peptide	3,93E-02	4,37
APOE_HUMAN	26	Apolipoprotein E	1,17E-03	3,61
APOC4_HUMAN	5	Apolipoprotein C-IV	2,27E-05	3,60
APOA4_HUMAN	28	Apolipoprotein A-IV	1,83E-05	3,54
CXCL7_HUMAN	6	Platelet basic protein	5,79E-03	3,05
PROS_HUMAN	13	Vitamin K-dependent protein S	2,07E-03	2,87
PLEK_HUMAN	4	Pleckstrin	1,34E-02	2,69
SEPP1_HUMAN	5	Selenoprotein P	3,31E-03	2,65
C4BPB_HUMAN	2	C4b-binding protein beta chain	4,79E-02	2,54
APOA1_HUMAN	22	Apolipoprotein A-I	2,38E-04	2,49
LBP_HUMAN	4	Lipopolysaccharide-binding protein	3,72E-02	2,38
ANT3_HUMAN	26	Antithrombin-III	6,32E-04	2,33
THRB_HUMAN	25	Prothrombin	9,76E-04	2,31
DESP_HUMAN	4	Desmoplakin	1,97E-02	2,29
FA12_HUMAN	12	Coagulation factor XII	1,15E-03	2,29
C4BPA_HUMAN	19	C4b-binding protein alpha chain	9,04E-03	2,25
FHR5_HUMAN	6	Complement factor H-related protein 5	7,01E-03	2,16
FA5_HUMAN	7	Coagulation factor V	2,80E-02	2,12
RAP1B_HUMAN	2	Ras-related protein Rap-1b	2,17E-02	2,00
CO4A_HUMAN	60	Complement C4-A	2,94E-02	1,95
CO4B_HUMAN	60	Complement C4-B	2,94E-02	1,95
GELS_HUMAN	19	Gelsolin	4,75E-04	1,94
PROP_HUMAN	3	Properdin	4,48E-02	1,87
IBP4_HUMAN	2	Insulin-like growth factor-binding protein 4	1,47E-02	1,84
APOC2_HUMAN	5	Apolipoprotein C-II	1,99E-02	1,83
MOES_HUMAN	2	Moesin	3,81E-02	1,80
PON1_HUMAN	10	Serum paraoxonase/arylesterase 1	3,35E-03	1,75
PPIB_HUMAN	2	Peptidyl-prolyl cis-trans isomerase B	2,51E-02	1,73
LIMS1_HUMAN	2	LIM and senescent cell antigen-like-containing domain protein 1	4,01E-02	1,73
CETP_HUMAN	2	Cholesteryl ester transfer protein	3,03E-02	1,53
QSOX1_HUMAN	4	Sulfhydryl oxidase 1	1,92E-02	1,50


Comparison GMI vs ST

Proteins with lower affinity with GMI implants





			GMI/ST	
			p value	Ratio
CO8B_HUMAN	10	Complement component C8 beta chain	3,64E-02	0,64
IGJ_HUMAN	3	Immunoglobulin J chain	4,59E-02	0,62
HV320_HUMAN	4	Ig heavy chain V-III region GAL	2,57E-02	0,62
TTHY_HUMAN	5	Transthyretin	1,70E-02	0,61
HBA_HUMAN	2	Hemoglobin subunit alpha	8,32E-03	0,53
ITIH4_HUMAN	18	Inter-alpha-trypsin inhibitor heavy chain H4	6,49E-03	0,51
CO8A_HUMAN	6	Complement component C8 alpha chain	2,48E-02	0,51
CO5_HUMAN	22	Complement C5	4,78E-02	0,49
IGHG2_HUMAN	14	Ig gamma-2 chain C region	1,38E-02	0,47
IGHG4_HUMAN	8	Ig gamma-4 chain C region	1,42E-02	0,44
FETUA_HUMAN	7	Alpha-2-HS-glycoprotein	1,08E-03	0,40
PGRP2_HUMAN	8	N-acetylmuramoyl-L-alanine amidase	1,81E-02	0,37
KNG1_HUMAN	17	Kininogen-1	6,98E-05	0,33
CFAB_HUMAN	19	Complement factor B	2,06E-02	0,33
FA11_HUMAN	20	Coagulation factor XI	6,27E-03	0,33
AMBP_HUMAN	8	Protein AMBP	5,48E-03	0,32
RET4_HUMAN	5	Retinol-binding protein	1,62E-03	0,31
VTDB_HUMAN	16	Vitamin D-binding protein	1,34E-02	0,30
IGHA1_HUMAN	11	Ig alpha-1 chain C region	1,34E-03	0,28
FCN1_HUMAN	2	Ficolin-1	7,85E-03	0,26
IGHA2_HUMAN	8	Ig alpha-2 chain C region	5,53E-03	0,24
HPT_HUMAN	21	Haptoglobin	4,07E-03	0,22
CERU_HUMAN	24	Ceruloplasmin	3,71E-04	0,22
CBG_HUMAN	2	Corticosteroid-binding globulin	5,97E-03	0,22
ANGT_HUMAN	6	Angiotensinogen	2,07E-03	0,21
AACT_HUMAN	11	Alpha-1-antichymotrypsin	4,84E-04	0,20
ITIH1_HUMAN	10	Inter-alpha-trypsin inhibitor heavy chain H1	2,50E-03	0,19
ALBU_HUMAN	57	Serum albumin	1,13E-02	0,18
A1AT_HUMAN	20	Alpha-1-antitrypsin	1,05E-03	0,17
TRFE_HUMAN	42	Serotransferrin	5,42E-03	0,15
A1AG2_HUMAN	7	Alpha-1-acid glycoprotein 2	2,30E-02	0,12
A1AG1_HUMAN	8	Alpha-1-acid glycoprotein 1	3,51E-03	0,12
ZA2G_HUMAN	11	Zinc-alpha-2-glycoprotein	7,13E-04	0,12
AFAM_HUMAN	10	Afamin	1,22E-03	0,12
A2MG_HUMAN	50	Alpha-2-macroglobulin	3,84E-03	0,11
HEMO_HUMAN	16	Hemopexin	4,91E-03	0,10
ITIH3_HUMAN	2	Inter-alpha-trypsin inhibitor heavy chain H3	1,47E-03	0,10
A1BG_HUMAN	9	Alpha-1B-glycoprotein	2,28E-03	0,09
A2GL_HUMAN	4	Leucine-rich alpha-2-glycoprotein	6,02E-03	0,08

Complex process





1º: IMMUNE RESPONSE





2º: BLOOD COAGULATION






Comparison **GMI** vs **ST**

GMI: Higher coagulation activity:

Pro- and anti- coagulant proteins were found more attached onto GMI surface




Accession	Description	p value	Ratio
PROC_HUMAN	Vitamin K-dependent protein C	1,17E-06	27,21
FA10_HUMAN	Coagulation factor X	1,86E-04	17,21
FA7_HUMAN	Coagulation factor VII	1,04E-04	16,27
FA9_HUMAN	Coagulation factor IX	2,44E-04	5,47
PROS_HUMAN	Vitamin K-dependent protein S	2,07E-03	2,87
ANT3_HUMAN	Antithrombin-III	6,32E-04	2,33
THRB_HUMAN	Prothrombin	9,76E-04	2,31
FA12_HUMAN	Coagulation factor XII	1,15E-03	2,29
FA5_HUMAN	Coagulation factor V	2,80E-02	2,12
KNG1_HUMAN	Kininogen-1	6,98E-05	0,33
FA11_HUMAN	Coagulation factor XI	6,27E-03	0,33
A2MG_HUMAN	Alpha-2-macroglobulin	3,84E-03	0,11
FA11_HUMAN	Coagulation factor XI	6,27E-03	0,33
A2MG_HUMAN	Alpha-2-macroglobulin	3,84E-03	0,11



PRO

REGULATION





ANTI

GMI implants IMPROVE coagulation process

Prothrombin


This protein display a key role in blood clotting




However, blood coagulation is a high regulated system.

ANT3, PROC, PROS could ensure the correct process development in order to avoid an acute inflammation and thrombosis

3º: OSTEOGENESIS



GMI IMPLANTS IMPROVE OSTEOGENESIS




Osteogenic proteins, higher affinity to GMI implants

GMI



Accession	Description
APOE_HUMAN	Apolipoprotein E
SEPP1_HUMAN	Selenoprotein P



GMI

Comparison GMI vs ST



Contents lists available at ScienceDirect

Biochimica et Biophysica Acta

journal homepage: www.elsevier.com/locate/bbagen



Review

Selenoproteins and selenium status in bone physiology and pathology



Zhichao Zhang^a, Jinsong Zhang^{b,*}, Jianru Xiao^{a,**}

^a Department of Orthopaedic Oncology, Changzheng Hospital, The Second Military Medical University, Shanghai 200003, PR China
^b School of Tea Food Science, Anhui Agricultural University, Hefei 230036, Anhui, PR China

ARTICLE INFO

Article history:
 Received 10 April 2014
 Received in revised form 23 July 2014
 Accepted 4 August 2014

ABSTRACT

Background: Emerging evidence supports the view that selenoproteins are essential for maintaining bone health. **Scope of review:** The current state of knowledge concerning selenoproteins and Se status in bone physiology and pathology is summarized.

Major conclusions: Antioxidant selenoproteins including glutathione peroxidase (GPx) and thioredoxin reductase (TrxR), as a whole, play a pivotal role in maintaining bone homeostasis and protecting against bone loss. GPx1, a major antioxidant enzyme in osteoclasts, is up-regulated by estrogen, an endogenous inhibitor of osteoclastogenesis. TrxR1 is an immediate early gene in response to 1 α ,25-dihydroxyvitamin D₃, an osteoblastic differentiation agent. The combination of 1 α ,25-dihydroxyvitamin D₃ and Se generates a synergistic elevation of TrxR1 activity in Se-deficient osteoblasts. Of particular concern, pleiotropic TrxR1 is implicated in promoting NF κ B activation. Coincidentally, TrxR inhibitors such as curcumin and gold compounds exhibit potent osteoclastogenesis inhibitory activity. Studies in patients with the mutations of selenocysteine insertion sequence-binding protein 2, a key trans-acting factor for the co-translational insertion of selenocysteine into selenoproteins have clearly established a causal link of selenoproteins in bone development. Se transport to bone relies on selenoprotein P. Plasma selenoprotein P concentrations have been found to be positively correlated with bone mineral density in elderly women.

General significance: A full understanding of the role and function of selenoproteins and Se status on bone physiology and pathology may lead to effectively prevent against or modify bone diseases by using Se.

Accession

APOE_HUMAN
 SEPP1_HUMAN

Description


Apolipoprotein E
 Selenoprotein P



Contents lists available at ScienceDirect

Bone

journal homepage: www.elsevier.com/locate/bone



Original Full Length Article

The role of Apolipoprotein E in bone metabolism

Andreas Niemeier^{a,b,*}, Thorsten Schinke^c, Joerg Heeren^b, Michael Amling^c

^a Department of Orthopaedics, University Medical Center Hamburg-Eppendorf, Hamburg, Germany

^b Department of Biochemistry and Molecular Biology II: Molecular Cell Biology, University Medical Center Hamburg-Eppendorf, Hamburg, Germany

^c Department of Osteology and Biomechanics, University Medical Center Hamburg-Eppendorf, Hamburg, Germany

ARTICLE INFO

Article history:
 Received 11 March 2011
 Revised 23 June 2011
 Accepted 8 July 2011
 Available online 23 July 2011

ABSTRACT

Apolipoprotein E (apoE) is a major structural apolipoprotein of several lipoprotein classes. Over the last 13 years, numerous studies have focused on the question whether human apoE affects bone phenotypes and, more recently, whether apoE regulates bone metabolism in mice. Here, we first provide a brief background introduction into the structure, established physiological and pathophysiological functions of apoE, and will then discuss the new aspects of the emerging role of apoE in bone.

Osteogenic potential of ApoE protein

Increase osteoblastic differentiation



Decrease osteoclastic differentiation



54

BioScience Trends. 2016; 10(1):54-66.

Original Article

DOI: 10.5582/bst.2016.01006

Multifarious effects of 17- β -estradiol on apolipoprotein E receptors gene expression during osteoblast differentiation *in vitro*

Yuyan Gui^{1,2,*}, Zhongliang Duan^{1,2,*}, Xuemin Qin^{1,2}, Wei Tang³, Hans-Jürgen Guber⁴, Dajin Li¹, Ling Wang^{1,2,**}

¹ Laboratory for Reproductive Immunology, Hospital & Institute of Obstetrics and Gynecology, IBS, Fudan University Shanghai Medical College, Shanghai, China;

² Shanghai Key Laboratory of Female Reproductive Endocrine Related Diseases, Shanghai, China;

³ Hepato-Biliary-Pancreatic Surgery Division, Department of Surgery, Graduate School of Medicine, The University of Tokyo, Tokyo, Japan;

⁴ Department of Pharmacy, Wagner Jauregg Hospital and Children's Hospital, Linz, Austria.

Summary

Apolipoprotein E (ApoE) regulated bone metabolism in mice might mediate uptake of lipid particles into target cells such as osteoblasts via receptor-mediated endocytosis by apoE receptors, which includes the low-density lipoprotein receptor (LDLR) family and heparan sulfate proteoglycans (HSPGs). There is no report regarding the expression of ApoE receptors mRNA induced by estrogen during osteoblast differentiation *in vitro*. Primary osteoblasts were collected from the calvaria of newborn mice and were subjected to osteoblast mineralization culture with serial concentrations of 17- β -estradiol (E2) *in vitro*. RNA was isolated at days 0, 5 and 25 of differentiation. Real-time PCR was conducted to analyze apoE receptors mRNA levels. We found that most LDLR family members genes were induced during osteoblast differentiation *in vitro*. The effect of E2 on apoE receptors gene expression during osteoblast differentiation was multifarious. The most noted members of the LDLR family involved in the maintenance of bone metabolism were LRP5, LRP6, LRP4, and Apoer2. LRP6 was up-regulated, while LRP5, LRP4, and Apoer2 were down-regulated by E2. Given that LRP6 is required for early stages of differentiation, we speculate E2 promotes osteoblast differentiation mainly in the early stage.


Keywords: 17- β -estradiol, Apolipoprotein E receptors, Low-density Lipoprotein Receptors Family, Heparan sulfate proteoglycans, Osteoblast differentiation, Reproductive endocrine metabolic network



Available online at www.sciencedirect.com

SciVerse ScienceDirect

journal homepage: www.elsevier.com/locate/yexcr



Research Article

Apolipoprotein E inhibits osteoclast differentiation via regulation of c-Fos, NFATc1 and NF- κ B

Woo-Shin Kim^a, Hyung Joon Kim^a, Zang Hee Lee^a, Youngkyun Lee^b, Hong-Hee Kim^{a,*}

^aDepartment of Cell and Developmental Biology, BK21 Program and Dental Research Institute, Seoul National University, 28 Yeongon-Dong, Chongno-Gu, Seoul 110-749, Republic of Korea

^bDepartment of Biochemistry, School of Dentistry, Kyungpook National University, Daegu 700-422, Republic of Korea

ARTICLE INFORMATION

Article Chronology:

Received 18 July 2012

Received in revised form

15 November 2012

Accepted 5 December 2012

Available online 12 December 2012

Keywords:

Osteoclast differentiation

Apolipoprotein E

c-Fos


NFATc1

NF- κ B

ABSTRACT

Apolipoprotein E (ApoE) plays a major role in the transport and metabolism of lipid. Other functions of ApoE include modulation of innate and adaptive immune responses. The expression of ApoE in osteoblasts and its relevance with bone formation have also been reported. However, the effect of ApoE on osteoclasts has not yet been examined. Here, we investigated the role of ApoE in osteoclast differentiation using bone marrow-derived macrophages (BMMs) and RAW264.7 cells. We found a down-regulation of ApoE gene expression during osteoclastic differentiation of those cells. Overexpression of ApoE in BMMs and RAW264.7 cells significantly blocked the induction of c-Fos and nuclear factor of activated T cell c1 (NFATc1), transcription factors critical for expression of osteoclast marker genes, by receptor activator of nuclear factor κ B ligand (RANKL), the osteoclast differentiation factor. ApoE inhibited osteoclast differentiation, as measured by decreased number of tartrate-resistant acid phosphatase (TRAP)-positive multinuclear cells (MNCs). In addition, ApoE reduced the expression of dendritic cell-specific transmembrane protein (DC-STAMP) and ATPase, H⁺ transporting, lysosomal 38 kDa, V0 subunit d2 (ATP6v0d2), genes involved in cell-cell fusion during osteoclastogenesis. Knock-down of ApoE using a specific siRNA promoted the RANKL-mediated induction of osteoclast differentiation. While ApoE did not affect the activation of ERK, JNK, and p38 MAPK signaling pathways by RANKL, the phosphorylation of p65 trans-activation domain on serine 536 and transcription activity of NF- κ B were reduced by ApoE overexpression. These findings suggest that ApoE plays an inhibitory role in osteoclast differentiation via the suppression of RANKL-dependent activation of NF- κ B and induction of c-Fos and NFATc1.

Role of vitamin K and ApoE: importance in bone metabolism



Rodrigues, A. *et al.* Low osteocalcin/collagen type I bone gene expression ratio is associated with hip fragility fractures. *Bone* 51, 981–9 (2012)

CONCLUSIONS

ACTIVITY	GMI	STRAUMMAN
IMMUNE RESPONSE	—	—
BLOOD COAGULATION	↑	↓
OSTEOGENESIS	↑	↓