SISTEMAS DE IMPLANTES IMPLANT SYSTEMS

www.gmidental.com

avantgard

avantgard
PEAK

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- Kit quirúrgico **gmi frontier/avantgard**
- Kit quirúrgico **gmi PEAK**

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KEY FEATURES

- Surface treatment / studies 4
- Recommendations of use 6
- Labelling description 7

IMPLANT SYSTEMS

- **gmi avantgard** dental implants 8
- Insertion sequence for **gmi avantgard** implants 9
- **gmi avantgard** implants range 10
- Drilling sequences for **gmi avantgard** implants 11
- **gmi avantgard PEAK** dental implants 13
- Insertion sequence for **gmi avantgard PEAK** implants 14
- **gmi avantgard PEAK** implants range 15
- Drilling sequences for **gmi avantgard PEAK** implants 16

ATTACHMENTS

- **gmi avantgard** healing abutments 17
- **gmi avantgard** implant carrier abutments 18
- **gmi avantgard** temporary cylinders 18
- **gmi avantgard** angled abutments 19
- **gmi avantgard** straight abutments 20
- **gmi avantgard** titanium base 21
- **gmi avantgard** multi-esthetic abutments 22
- **gmi avantgard** bioesthetic bases 25
- **gmi avantgard** CrCo pre-milled blanks 26
- **gmi avantgard** replicas 27
- **gmi avantgard** impression coping 28
- **gmi avantgard** clinic and laboratory screws 29
- **gmi avantgard** MEA and BE screws 30
- Recommended tightening torques 31

GMI SURGICAL KITS

- **gmi frontier/avantgard** surgical kit 33
- **gmi PEAK** surgical kit 34

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


To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.





SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

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The **gmi avantgard** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

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The **gmi avantgard** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi avantgard** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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


To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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
In the chart are detailed all the symbols appearing on the gmi avantgard implant system's labelling and packaging and their corresponding description.

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* Minimum requirements for downloading and opening the electronic instructions for use:
- Stable internet connectivity.
- A device with a modern operating system that supports PDF viewing.
- Adequate storage space on your device.
- A web browser such as Google Chrome, Mozilla Firefox, Safari, Microsoft Edge or Opera installed on your device.
- A PDF reader application installed on your device..

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
gmi avantgard conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high primary stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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gmi avantgard conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi avantgard** implants consists of a titanium support which houses the implant-implant carrier set, prepared to be inserted directly into the mouth, using a single hex-3.00 mm implant carrier wrench, designed for this purpose. This new system offers the dentist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.*

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* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the carton.
2. **Remove the secondary packaging and product identification labels.**
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. **Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.**
5. Extract primary plug making a lateral movement, and save it for posterior use.
6. **Insert the implant carrier wrench until a slight retention is noticed, addressing the notches on the key with the faces of the implant carrier hex broaching.**
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. **Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the subcrestal level. Remove the key.**
- 9-10. Hold the implant carrier with flat manual wrench and untighten screw using HEX-1.20 mm wrench until the cone is disassembled. Remove the implant carrier.
11. **Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.**
12. **Manual tightening cover screw to the implant. Remove key and suture the incision.**

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


The following is a description of the drilling sequences for the different conical connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped drill: 1200 - 1500 rpm.
- Pilot drill: 700 - 900 rpm.
- Step drills:
 - ∅3.00 & ∅3.50 mm → 400 - 700 rpm.
 - ∅4.00 & ∅5.10 mm → 400 - 600 rpm.
- Cortical drills: 200 - 400 rpm.

For drilling, external cooling with saline solution is required.


avantgard ø3,75			
Descripción - Description		Ref.	
1	Fresa lanza	Lance-shaped drill	KYFOC1225
2	Fresa piloto Ø2,50	Pilot drill Ø2.50	KYFOC2250
3	Fresa escalonada Ø2,50/Ø3,00	Step drill Ø2.50/Ø3.00	KYFOC5130
4	Fresa cortical Ø4,00	Cortical drill Ø4.00	KYFOC4340
5	Fresa escalonada Ø3,00/Ø3,50	Step drill Ø3.00/Ø3.50	KYFOC5135


avantgard ø4,25			
Descripción - Description		Ref.	
1	Fresa lanza	Lance-shaped drill	KYFOC1225
2	Fresa piloto Ø2,50	Pilot drill Ø2.50	KYFOC2250
3	Fresa escalonada Ø2,50/Ø3,00	Step drill Ø2.50/Ø3.00	KYFOC5130
4	Fresa escalonada Ø3,00/Ø3,50	Step drill Ø3.00/Ø3.50	KYFOC5135
5	Fresa cortical Ø4,50	Cortical drill Ø4.50	KYFOC4345
6	Fresa escalonada Ø3,50/Ø4,00	Step drill Ø3.50/Ø4.00	KYFOC5140

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
gmi avantgard PEAK conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the conical core and progressive thread profile facilitate implant insertion and act as a bone compactor allowing high primary stability to be obtained in sparse bones.

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gmi avantgard PEAK conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi avantgard PEAK** implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single Avantgard RP implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.*

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* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the RP insertion wrench facing the grooves of the wrench with those of the implant. Apply moderate pressure to drive the wrench into the implant well.
7. Check that the wrench is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the subcrestal level. Remove the wrench with gentle lateral movements.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove wrench and suture the incision.

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



Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

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Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Available in emergency heights from 1.50 to 4.50 mm.
- Apical bit to facilitate the start of threading.
- Manual tightening torque: Max. 15 N-cm.

Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
RP ●●	Ø 5,00	1,50	KDPoC4001
		2,50	KDPoC4002
		3,50	KDPoC4003
		4,50	KDPoC4004

Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
RP ●●	Ø 6,00	2,50	KDPoC4009
		3,50	KDPoC4010
		4,50	KDPoC4011


Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
WP ●	Ø 6,00	1,50	KDPoC4005
		2,50	KDPoC4006
		3,50	KDPoC4007
		4,50	KDPoC4008

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When the use of a custom abutment is required, implant carrier abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, implant carrier abutments can be used like impression coping, threading them to the implant screw through the short impression screw and using them to make the negative of the patient's mouth for subsequent positation on plaster.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 3.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Premounted clinic screw included.
- Screw tightening torque: 25 N-cm.


Plataforma Platform	Ø hombro Ø shoulder	Ref.
RP ●●	Ø 3,75 - Ø 4,25	KDBoF4001
WP ●	Ø 4,75	KDBoF4002
	Ø 5,75	KDBoF4003

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When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary abutment is used, defined as attachment that once attached to the crown of the implant with clinic screw, will absorb most of the masticatory stresses while avoiding overloading the implant while the osseointegration process is completed.

Key features

- Machined in high strength technical plastic.
- Rotary version only for multiple restorations.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw included.
- Manual tightening torque: Max. 15 N-cm.

Plataforma Platform	Tipo conexión Connection type	Ref.
RP ●●	Anti-rotatoria • Non-rotary	KDHoF4001
	Rotatoria • Rotary	KDHoF4002
WP ●	Anti-rotatoria • Non-rotary	KDHoF4003
	Rotatoria • Rotary	KDHoF4004

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

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Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw included.
- Abutment shoulder adapted to the shape of the gum.
- Screw tightening torque: 25 N-cm.

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP ●●	15°	2,50	KDMoF4013
	15°	4,50	KDMoF4014
	20°	2,50	KDMoF4015
	20°	4,50	KDMoF4016

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
WP ●	15°	2,50	KDMoF4017
	15°	4,50	KDMoF4018
	20°	2,50	KDMoF4019
	20°	4,50	KDMoF4020

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

Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Available only in non-rotary version.
- Clinic screw included.
- Available in four emergency heights from 1.50 to 4.50 mm.
- Screw tightening torque: 25 N-cm.

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
RP ●●	Ø 4,00	1,50	KDMoF4001
		2,50	KDMoF4002
		3,50	KDMoF4003
		4,50	KDMoF4004
	Ø 5,00	1,50	KDMoF4005
		2,50	KDMoF4006
		3,50	KDMoF4007
		4,50	KDMoF4008

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
WP ●	Ø 6,00	1,50	KDMoF4009
		2,50	KDMoF4010
		3,50	KDMoF4011
		4,50	KDMoF4012

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
Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw included.
- Available in emergency height of 1.50 mm and 3.00 mm.
- Non-rotary titanium base with indexation.
- Screw tightening torque: 25 N-cm.
- Scan directly on Titanium base or alternative scan on Scan abutment.

Base titanio - Titanium base

Plataforma Platform	Altura (mm) Height (mm)	Tipo conexión Connection type	Ref.
RP ● ●	1,50	Anti-rotatoria • Non-rotary	KDMoF4021
	1,50	Rotatoria • Rotary	KDMoF4023
	3,00	Anti-rotatoria • Non-rotary	KDMoF4025
	3,00	Rotatoria • Rotary	KDMoF4027
WP ●	1,50	Anti-rotatoria • Non-rotary	KDMoF4022
	1,50	Rotatoria • Rotary	KDMoF4024
	3,00	Anti-rotatoria • Non-rotary	KDMoF4026
	3,00	Rotatoria • Rotary	KDMoF4028

Calcínable - Castable

Plataforma Platform	Altura total (mm) Total Height (mm)	Cantidad (Ud.) Quantity (Un.)	Ref.
RP ● ●	11,00	10	KDCoF3605
WP ●	11,00	10	KDCoF3606

Base titanio gmi avantgard - gmi avantgard titanium base

Pilar escaneado - Scan abutment

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP ●●	6,30	KDioF4001
WP ●	6,30	KDioF4002

Cuerpo pilar escaneado - Scanbody

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP - WP ●●	7,50	KDioC4003

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.50 to 4.50 mm.
- One-piece straight abutments in rotary version.
- Tightened with hex wrench: B/F 2.00 mm (KYLOCo149).
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.
- Tightening torques: ME straight abut. (30 N-cm), ME angled abut. (20 N-cm) and healing caps (Manual: Max. 15 N-cm).

Pilares multi-estético gmi avantgard - gmi avantgard multi-esthetic abutments

Pilares multi-estético inclinados Multi-esthetic angled abutments

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
	17°	2,50	KDGoF4009
RP ●●	17°	3,50	KDGoF4010
	30°	4,00	KDGoF4011

Pilares multi-estético rectos Multi-esthetic straight abutments

Plataforma Platform	Altura (mm) Height (mm)	Ref.
	1,50	KDGoF4001
RP ●●	2,50	KDGoF4002
	3,50	KDGoF4003
	4,50	KDGoF4004
WP ●	1,50	KDGoF4005
	2,50	KDGoF4006
	3,50	KDGoF4007
	4,50	KDGoF4008

Pilares cicatrización PME Healing abutments MEA


Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP - WP ●●	Ø 4,80	KDPOF0035
	Ø 6,00	KDPOF0037

Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

Plataforma Platform	Ref.
RP - WP ●●	KDS0F0015

Pilares multi-estético gmi avantgard - gmi avantgard multi-esthetic abutments

Tomas impresión cubeta cerrada (CC) PME
 Closed tray (CT) impression coping MEA

Plataforma Platform	Cantidad (Ud.) Quantity (Un.)	Ref.
RP - WP ●●	1	KDS0F0017
	5	KDS0F5001

Cilindro provisional PME
 Temporary cylinder MEA


Plataforma Platform	Ref.
RP - WP ●●	KDGoF0028

Calcinable PME
 Castable MEA

Plataforma Platform	Ref.
RP - WP ●●	KDCoF0117

Base titanio PME
 Titanium base MEA

Plataforma Platform	Ref.
RP - WP ●●	KDGoF0035

Scanbody de titanio PME
 Titanium scanbody MEA

Plataforma Platform	Ref.
RP - WP ●●	KDIOF0002

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Once the insertion of the implant in the bone bed has been completed, the bioesthetic base is used, which is defined as the machined attachment that, once fixed to the coronal part of the implant, allows single restorations at the abutment level.

Placing the bioesthetic base during surgery allows the soft tissues to heal around it.

Once the tissue has healed, it serves as the basis for the rehabilitation, avoiding the connection and disconnection of the components to the implant and, therefore, keeping the biological seal intact.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 2.50 to 4.50 mm.
- Straight two-piece bases, including the screw.
- Tightened with hex wrench: B/F 2.00 mm (KYLOCo149).
- Impression coping with screw included.
- Castable.
- Tightening torques: BE base (25 N-cm) and healing abutment (Manual: Max. 15 N-cm).

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
CrCo pre-milled blanks consist of a block with the implant connection machined, which allows an excellent fit between the implant and the abutment.

They are indicated to produce custom one-piece abutments using CAD/CAM allowing direct application of ceramic.

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

Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

Key features


- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.


Réplicas de implante Implant replicas	
Plataforma Platform	Ref.
RP ●●	KDRoC4001
WP ●	KDRoC4002

Réplicas digitales Digital replicas	
Plataforma Platform	Ref.
RP ●●	KDRoF4003
WP ●	KDRoF4004

Réplicas para pilares multi-estético MEA replicas	
Plataforma Platform	Ref.
RP - WP ●●	KDRoC0012

Réplicas digitales para pilares multi-estético MEA digital replicas	
Plataforma Platform	Ref.
RP - WP ●●	KDRoF0014

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
Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, is used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- Manual tightening torque: Max. 15 N-cm.

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
Clinic screws are defined as machined devices used for the final anchoring of the abutments or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of abutments to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

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Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo clínica Clinic screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP ●●	⊕ 1,20 mm	KDToC4003
	⊗ TX	KDToC4021
WP ●	⊕ 1,20 mm	KDToC4004
	⊗ TX	KDToC4022
Tornillo laboratorio Laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP ●●	⊕ 1,20 mm	KDToC4005
	⊗ TX	KDToC4023
WP ●	⊕ 1,20 mm	KDToC4006
	⊗ TX	KDToC4024

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MEA and BE clinical screws are defined as machined devices used for the final anchoring of the ME and BE abutments, while laboratory screws are the machined devices used for the temporary anchoring of ME and BE abutments to the replica during the fabrication of the prostheses in the laboratory.

Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm and 2.00 mm hex wrench.
- TX screw allows the screw channel to be angled up to 20°.

Tornillo PME MEA screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ● ●	○ 1,20 mm	KDToC0056
	⊛ TX	KDToC0070

Tornillo laboratorio PME MEA laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ● ●	○ 1,20 mm	KDToC0057
	⊛ TX	KDToC0072

Tornillo BE BE screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ● ●	○ 1,20 mm	KDToC4016
	⊛ TX	KDToC0074

Tornillo laboratorio BE BE laboratory screw		
Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP - WP ● ●	○ 1,20 mm	KDToC4020
	⊛ TX	KDToC0075

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

gmi frontier/avantgard surgical kit consist of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi frontier** and **avantgard** dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134 °C.

To complement this kit, **gmi** has the **PEAK drill Kit**, allowing it to be easily expanded for the placement of 2 additional implant systems (**frontier PEAK** and **avantgard PEAK**).

Kit quirúrgico gmi frontier/avantgard - gmi frontier/avantgard surgical kit

Descripción - Description		Un.	Ref.
Fresa lanza	Lance-shaped drill	1	KYFoC1225
Fresa bola ø2,0 mm	Ball drill ø2.0 mm	1	KYFoC1502
Fresa piloto ø2,50 mm	Pilot drill ø2.50 mm	1	KYFoC2250
Fresa escalonada ø2,50 / ø3,00 mm	Step drill ø2.50 / ø3.00 mm	1	KYFoC5130
Fresa escalonada ø3,00 / ø3,50 mm	Step drill ø3.00 / ø3.50 mm	1	KYFoC5135
Fresa escalonada ø3,50 / ø4,00 mm	Step drill ø3.50 / ø4.00 mm	1	KYFoC5140
Fresa escalonada ø4,00 / ø4,50 mm	Step drill ø4.00 / ø4.50 mm	1	KYFoC5145
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4336
Fresa cortical ø4,00 mm	Cortical drill ø4.00 mm	1	KYFoC4340
Fresa cortical ø4,50 mm	Cortical drill ø4.50 mm	1	KYFoC4345
Fresa cortical ø4,90 mm	Cortical drill ø4.90 mm	1	KYFoC4349
Fresa cortical ø5,90 mm	Cortical drill ø5.90 mm	1	KYFoC4359
Llave transportador plana 45°	Manual wrench 45°	1	KYLoC0093
Llave carraca IP	TI ratchet wrench	1	KYLoF0113
Llave-transportador HEX-2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYLoF0096
Llave-transportador HEX-2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYLoF0111
Llave-transportador HEX-2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYLoF0091
Llave-transportador Avantgard HEX-3,00 mm IP corta	Short HEX-3.00 mm Avantgard implant carrier TI wrench	1	KYLoF0137
Llave-transportador Avantgard HEX-3,00 mm IP larga	Long HEX-3.00 mm Avantgard implant carrier TI wrench	1	KYLoF0138
Llave-transportador Avantgard HEX-3,00 mm CA	HEX-3.00 mm Avantgard implant carrier HP wrench	1	KYLoF0136
Llave carraca corta HEX-1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLoF0128
Llave carraca larga HEX-1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLoF0129
Alargador fresas	Drill extender	1	KYLoC0023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLoC0078
Indicador profundidad ø2,50 mm	Depth gauge ø2.50 mm	1	KYLoC0165

Kit fresas gmi PEAK - gmi PEAK drill kit

Descripción - Description.		Un.	Ref.
Fresa escalonada ø2,00 / ø2,80 mm	Step drill ø2.00 / ø2.80 mm	1	KYFoC5128
Fresa escalonada ø3,20 / ø3,65 mm	Step drill ø3.20 / 3.65 mm	1	KYFoC5137
Fresa escalonada ø3,70 / ø4,15 mm	Step drill ø3.70 / ø4.15 mm	1	KYFoC5142
Fresa escalonada ø4,15 / ø4,65 mm	Step drill ø4.15 / ø4.65 mm	1	KYFoC5147
Fresa cónica ø3,20 mm	Conical drill ø3.20 mm	1	KYFoC7032
Fresa cónica ø3,70 mm	Conical drill ø3.70 mm	1	KYFoC7037
Fresa cónica ø4,20 mm	Conical drill ø4.20 mm	1	KYFoC7042
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFoC4337
Fresa cortical ø4,15 mm	Cortical drill ø4.15 mm	1	KYFoC4342
Fresa cortical ø4,65 mm	Cortical drill ø4.65 mm	1	KYFoC4347

Llaves de inserción avantgard PEAK vendidas por separado - Avantgard PEAK insertion wrenches sold separately

Kit quirúrgico gmi PEAK - gmi PEAK surgical kit

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gmi PEAK surgical kit consist of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi frontier PEAK** and **avantgard PEAK** dental implant systems. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi surgical kits are distributed in boxes made of high strength technical polymer, and must be sterilized in a moist heat autoclave before use at a maximum temperature of 134°C.

To complement this kit, **gmi** has the **frontier/avantgard drill Kit**, allowing it to be easily expanded for the placement of 2 additional implant systems (frontier and avantgard).

Kit quirúrgico gmi PEAK - gmi PEAK surgical kit

Descripción - Description		Un.	Ref.
Fresa bola ø2,0 mm	Ball drill ø2.0 mm	1	KYFOC1502
Fresa inicial ø2,0 mm	Initial drill ø2.0 mm	1	KYFOC1401
Fresa cónica ø3,20 mm	Conical drill ø3.20 mm	1	KYFOC7032
Fresa cónica ø3,70 mm	Conical drill ø3.70 mm	1	KYFOC7037
Fresa cónica ø4,20 mm	Conical drill ø4.20 mm	1	KYFOC7042
Fresa cortical ø3,65 mm	Cortical drill ø3.65 mm	1	KYFOC4337
Fresa cortical ø4,15 mm	Cortical drill ø4.15 mm	1	KYFOC4342
Fresa cortical ø4,65 mm	Cortical drill ø4.65 mm	1	KYFOC4347
Fresa escalonada ø2,00 / ø2,80 mm	Step drill ø2.00 / ø2.80 mm	1	KYFOC5128
Fresa escalonada ø3,20 / ø3,65 mm	Step drill ø3.20 / ø3.65 mm	1	KYFOC5137
Fresa escalonada ø3,70 / ø4,15 mm	Step drill ø3.70 / ø4.15 mm	1	KYFOC5142
Fresa escalonada ø4,15 / ø4,65 mm	Step drill ø4.15 / ø4.65 mm	1	KYFOC5147
Indicador profundidad ø2,00 mm	Depth gauge ø2.00 mm	1	KYLOCO164
Indicador profundidad cónico ø3,20 mm	Conical depth gauge ø3.20 mm	1	KYLOCO159
Indicador profundidad cónico ø3,70 mm	Conical depth gauge ø3.70 mm	1	KYLOCO160
Indicador profundidad cónico ø4,20 mm	Conical depth gauge ø4.20 mm	1	KYLOCO161
Llave carraca IP	TI ratchet wrench	1	KYLOFO113
Llave Avantgard CA corta	Short Avantgard HP wrench	1	KYLOCO151
Llave Avantgard carraca corta	Short Avantgard ratchet wrench	1	KYLOCO157
Llave Avantgard carraca larga	Long Avantgard ratchet wrench	1	KYLOCO162
Llave carraca corta HEX - 1,20 mm	Short HEX-1.20 mm ratchet wrench	1	KYLOFO128
Llave carraca larga HEX - 1,20 mm	Long HEX-1.20 mm ratchet wrench	1	KYLOFO129
Alargador fresas	Drill extender	1	KYLOCO023
Paralelizador ø2,00 / ø2,80 mm	Paralleling pin ø2.00 / ø2.80 mm	3	KYLOCO078
Llave-transportador HEX - 2,45 mm IP corta	Short HEX-2.45 mm implant carrier TI wrench	1	KYLOFO096
Llave-transportador HEX - 2,45 mm IP larga	Long HEX-2.45 mm implant carrier TI wrench	1	KYLOFO111
Llave-transportador HEX - 2,45 mm CA	HEX-2.45 mm implant carrier HP wrench	1	KYLOFO091

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